


IMPLANT LINE ISY KONE



GUIDELINES FOR THE USE
OF ISY KONE
IMPLANT SYSTEM

IMPORTANT NOTE

For latest updates and information, visit **www.btk.dental**

This manual provides dental practitioners and related specialists with general information regarding the use of ISY KONE dental implant systems.

For detailed information on other specific implant lines and their restorative procedures, please refer to the corresponding manuals, specific literature or refer to the BTK website.

Consider to regularly visit practical courses for updates and professional exchange with dedicated colleagues in order to ensure your long-term success with implant-borne dental restorations.

© 2021 BTK - the smile system.

ISY KONE

ONE BONE-LEVEL IMPLANT:
ALL PROSTHETIC OPTIONS

INDEX

1.	CORPORATE BACKGROUND	6
2.	IMPLANT SELECTION	8
3.	IMPLANT CHARACTERISTICS	10
4.	IMPLANT PORTFOLIO	12
5.	IMPLANT INSERTION PROCEDURE	14
6.	SURGICAL KIT, TORQUE WRENCH AND DRILLS	16
7.	SURGICAL PROTOCOLS	22
8.	ADDITIONAL INSTRUMENTS	24
9.	MORSE TAPER MTH KN/KR	26
10.	INTERNAL HEXAGON INT IR/IW	30
11.	EXTERNAL HEXAGON EXT EN/ER	35
12.	OVERDENTURE	40
13.	MATERIAL SPECIFICATIONS	42
14.	LABELS SYMBOLS	44
15.	DELIVERY TERMS AND CONDITIONS	45



CORPORATE BACKGROUND

Privately held BTK BIOTEC was founded in 1998 in order to improve the quality of life of people with missing teeth.

BTK is a dedicated supporter of the genuine “100% Made in Italy” label, because with this it is guaranteed that BTK products are of unmatchable Italian craftsmanship and premium quality materials offering dedicated specialization and ample differentiation.



BTK Headquarters - NORTH ITALY

Implanting Trust, Smile again!

By combining cutting-edge technologies and biology, BTK's mission is to offer affordable and personalized implant-borne solutions thereby sustainably improving the daily life of dental patients.

Together with leading professionals, BTK strives to become a reference in replacing missing teeth with trusted implant solutions in order to improve oral health around the globe.





PREMIUM QUALITY MATERIALS

Grade 4 commercially pure titanium (ASTM F 67 / ISO 5832-2) is BTK's material of choice for dental implants. Grade 4 is slightly harder to work, but it provides the highest strength and durability characteristics among the commercially pure titanium grades, making it the natural choice for BTK dental implants.

Grade 5 titanium (ASTM F 136 / ISO 5832-3) is used for BTK's prosthetic components, as these are subject to certain levels of stress and in the MINI line implants. This high-strength version, also known as Ti-6Al-4V, is widely used in orthopedics and shows excellent long-term physical and mechanical properties.



ENDOSSEOUS SURFACE DAE

Clinical trials confirm that roughened endosseous surfaces perform better than machined surfaces concerning endosseous wound healing, "de novo" bone formation and reduced time-to-loading.

Our DAE (dual- acid-etched) process aims to obtain a moderately rough surface with a controlled micro-roughness.



IMPLANT-ABUTMENT CONNECTION

The precision of the connection between implant and abutment creating a tight seal may be beneficial in preventing inflammatory bacteria propagating in the interface between different components.

Apart from that, extremely tight tolerances as applied by BTK help to avoid micromovements.

Providing precision in every part produced is one of our key contributions ensuring longterm restorative success.



RESTORATIVE OPTIONS

The purpose of dental implant therapy, now widely used in dentistry, is to replace lost dental elements with biocompatible titanium implants, in order to obtain a new and correct occlusion, using prostheses on implants.

In order to achieve this goal, BTK offers a focused portfolio of restorative solutions backed-up by comprehensive clinical experience. BTK offers a variety of prostheses components to satisfy the clinical preferences and needs of the patients.



MADE IN ITALY, USED GLOBALLY

We constantly ensure that the quality of our products and services meet the high expectations of our customers and their patients. Specialized professionals are taking care to offer comprehensive solutions in applied research, engineering, education and related activities.

Our brand is a solid promise of quality, we are certified UNI EN ISO 9001, UNI EN ISO 13485 and MDD 93/42/EEC and subsequent amendments and additions, and is therefore authorized to apply the CE Mark on its products.

IMPLANT PORTFOLIO


BTK is dedicated to offer comprehensive implant solutions to meet the requirements of individual clinical situations, user preferences and economic constraints.

Different designs, sizes, diameters, surfaces and abutment connections are available, while at the same time BTK strives to maintain a small number of precision-instruments thus simplifying procedures and limiting investments needed.



SELECTION OF THE IMPLANTS ON THE BASIS OF BONE DENSITY AND OF THE DRILLING PROTOCOL

IMPLANT DESIGN



CHARACTERIZATION OF BTK IMPLANT-ABUTMENT CONNECTIONS




MORSE-TAPER (MTH/MTD)

BTK's morse-taper hexagon connection comprises a 2.6 mm conical portion at 11° above a hexagon configuration combined with a M1.4 (KN), M1.6 (KR/DR) or M1.8 (KW) abutment screw to deliver adequate pre-load with a minimum of tightening.

Implants with a tapered interface can resist larger axial and transversal forces than implants with a flat interface. The design guides the abutment into a predictable location with a precise fit with the inner portion of the implant.

The precision of the conical connection with its tight seal may be beneficial in preventing inflammatory bacteria from propagating in the interface between implant and abutment and it helps to avoid micro-movements.




INTERNAL HEXAGON (INT)

BTK's internal hexagon connection comprises a parallel hexagon of 2 mm length opening with a small conical portion combined with a M1.8 abutment screw to deliver adequate pre-load with a minimum of tightening.

The internal hexagon has two functions: to transfer the torque momentum during implant placement and as an indexing system to transfer the precise 3D-position of the implant to the master cast.

Internal indexing systems have some advantages over external indexing systems since they allow longer engaging surfaces while reducing the platform height of the implant. This offers somewhat more flexibility in designing the emergence profile of the final restoration.



EXTERNAL HEXAGON (EXT)

BTK's external hexagon connection comprises a parallel hexagon at 0.7 mm height and a 90° shoulder to allow a flat-to-flat margin fit to the implant. Abutments are connected to the implant using a M1.8 (EN) or M2.0 (ER/EW) abutment screw.


The abutment screw plays a central role for the mechanical, long-term strength and fatigue resistance of the implant abutment connection. The requirements for such a screw are many, such as no loosening, long-term fatigue resistance, overload protection and safe pick-up and handling ability.

Due to the fact that the abutment screw is exposed to heavy dynamic loads, the precise application of tightening torque force is essential.

MORSE-TAPER (MTH)	INTERNAL HEXAGON (INT)	EXTERNAL HEXAGON (EXT)
KN = KONIC NARROW		EN = EXTERNAL NARROW
KR = KONIC REGULAR DR = DODECAGON REGULAR	IR = INTERNAL REGULAR	ER = EXTERNAL REGULAR
	IM = INTERNAL MEDIUM	
KW = KONIC WIDE	IW = INTERNAL WIDE	EW = EXTERNAL WIDE

NOTE that different BTK implants require different types of prosthetic platforms using corresponding abbreviations according to their sizes. For more details, refer to the corresponding BTK implant lines documentation.

3. IMPLANT CHARACTERISTICS




ONE ENDOSSEOUS DESIGN: ALL RESTORATIVE OPTIONS

ISY KONE implants are indicated for single or multiple unit restorations in splinted or non-splinted applications.

They are suitable for two-stage or single-stage healing combined with immediate, early or delayed loading protocols if sufficient primary stability and appropriate loading is ensured.

Restorations range from removable / fixed single crowns to partial / full dental arch restorations or overdenture applications to restore the masticatory function and appropriate aesthetics.

The key feature of the ISY KONE concept comprises a standardized endosseous screw-type implant body combined with three different prosthetic connections offering maximum versatility based on a single implant thread type and a corresponding surgical kit.





IMPLANT CHARACTERISTICS


This bone-level, root-shaped and self-tapping implant offers a back-tapered, micro-threaded smooth neck portion of 0,8-1,15 mm and is well suitable for most indications. ISY KONE is placed at bone level and offers both, healing modes. The single-lead square thread pitch on ISY KONE measures 0,9 mm for all diameters including Ø 2.9 mm / Ø 3.3 mm / Ø 3.7 mm / Ø 4.1 mm / Ø 4.8 mm / Ø 6 mm.

ISY KONE is a highly versatile type of implant line as it offers three different abutment connections either external (EXT) or internal (INT) hexagon, or as a third option the well-proven 11° (MTH) hexagon having regular (KR) and narrow (KN) corresponding prosthetic components. It shares the surgical kit with BT SAFE, BT NANO and ISY+.


PORTFOLIO ISY KONE

MORSE-TAPER CONNECTION KR - KN	implant length in mm B							APEX TIP Ø / mm D
	Ø mm A	8	10	12	14	16	18	
	KONIC REGULAR (KR) C	KR (konic regular) M1.6 occlusal thread / prosthetic platform Ø 2.8 mm / smooth neck portion 0.95 mm						
	3,3	129KR33J	129KR33L	129KR33N	129KR33Q	-	-	2.0
	3,7	129KR37J	129KR37L	129KR37N	129KR37Q	129KR37S	129KR37T	2.4
	4,1	129KR41J	129KR41L	129KR41N	129KR41Q	129KR41S	129KR41T	2.7
	4,8	129KR48J	129KR48L	129KR48N	129KR48Q	-	-	3.3
	6	129KR60J	129KR60L	-	-	-	-	4.3
	KONIC NARROW (KN) C	KN (konic narrow) M1.6 occlusal thread / prosthetic platform Ø 2.8 mm / smooth neck portion 0.5 mm						
	2,9	129KN29J	129KN29L	129KN29N	129KN29Q	-	-	1.7


INTERNAL HEXAGON CONNECTION IR / IW	implant length in mm B							APEX TIP Ø / mm D
	Ø mm A	8	10	12	14	16		
	INTERNAL HEXAGON (IR) C	IR (internal regular) M1.8 occlusal thread / prosthetic platform Ø 3.5 mm / smooth neck portion 0.95 mm						
	3,3	129IR33J	129IR33L	129IR33N	129IR33Q	-	-	2.0
	3,7	129IR37J	129IR37L	129IR37N	129IR37Q	129IR37S	-	2.4
	4,1	129IR41J	129IR41L	129IR41N	129IR41Q	129IR41S	-	2.7
	4,8	129IR48J	129IR48L	129IR48N	129IR48Q	-	-	3.3
	INTERNAL HEXAGON (IW) C	IW (internal wide) M1.8 occlusal thread / prosthetic platform Ø 5.0 mm / smooth neck portion 0.95 mm						
	4,8	129IW48J	129IW48L	129IW48N	129IW48Q	-	-	3.3

EXTERNAL HEXAGON CONNECTION EN / ER	implant length in mm B							APEX TIP Ø / mm D
	Ø mm A	8	10	12	14	16	18	
	EXTERNAL HEXAGON (EN) C	EN (external narrow) M1.8 occlusal thread / prosthetic platform Ø 3.4 mm / smooth neck portion 0.8 mm						
	3,3	131EN33J	131EN33L	131EN33N	131EN33Q	-	-	2.0
	EXTERNAL HEXAGON (ER) C	ER (external regular) M2.0 occlusal thread / prosthetic platform Ø 4.1 mm / smooth neck portion 1.15 mm (Ø 3.7) 0.8 mm (Ø 4.1 / Ø 4.8)						
	3,7	131ER37J	131ER37L	131ER37N	131ER37Q	131ER37S	131ER37T	2.4
	4,1	131ER41J	131ER41L	131ER41N	131ER41Q	131ER41S	131ER41T	2.7
	4,8	131ER48J	131ER48L	131ER48N	131ER48Q	-	-	3.3


**MORSE-TAPER
CONNECTION KR/KN**



**INTERNAL HEXAGON
CONNECTION IR / IW**



**EXTERNAL HEXAGON
CONNECTION EN / ER**



The color codes applied for different implant diameters and prosthetic platforms are indicated below:

	RED	PURPLE	WHITE	BLUE	YELLOW	GREEN
IMPLANT DIAMETER Ø	2,9	3,3	3,7	4,1	4,8	6
PROSTHETIC PLATFORM MTH	KN	KR	KR	KR	KR	KR
PROSTHETIC PLATFORM INT		IR	IR	IR	IR/IW	
PROSTHETIC PLATFORM EXT		EN	ER	ER	ER	

HANDLING OF STERILE IMPLANT PACKAGING

CAUTION

The sealed package of the medical device (MD) must be opened in a surgically suitable environment.

The removal of the implant and of the cover screw, if provided, must be carried out using sterilized instruments, avoiding any contact with non-sterile surfaces.

The sterility of the medical device is only guaranteed if the following conditions are met:

the expiry date stated on the packaging is still valid; there is a red dot on the sterile vial that signals the successful operation of gamma ray irradiation; the sealed package has not been opened and does not show damage or perforations. If only one of the aforementioned conditions is not respected, the device must not be used.

The device is disposable; the reuse can compromise the safety features of the device making it inappropriate for its intended use. BIOTEC explicitly declares that the MD is for single use and assumes no responsibility for any re-use by users.



1

BTK dental implants are supplied sterile in a double-vial package. The implant diameter, length and lot are shown on the label located in the vial containing the implant and in the outer label on the back of the packaging.



2

Open the blister from the back by breaking the outer label, and take out the vial.



3

The top lid of the vial is protected by the seal label. The color of the seal label identifies the diameter of the implant. To facilitate compliance with the traceability requirement of the medical device, there are two detachable patient-labels in the vial. One must be stuck onto the patient's medical record and one onto the patient's implant passport.



4

Open the external vial and withdraw the internal vial containing the implant in a surgically suitable environment. The internal vial must be handled with sterile gloves.



5

Remove the safety cap of the sterile inner vial, which always includes the sterile closure screw. **WARNING** The internal vial consists of 3 parts. The cover screw (locking screw), if provided, is placed in the vial cap. Hold the vial upright to prevent the devices from leaking out. Unscrew the central part of the vial, to access the implant.



6

Some implant lines are supplied with mounting device connected to the implant, other lines are supplied without. Depending on the different configuration, use the appropriate instrument for the implant withdrawal from the vial and for the insertion of the same in the previously prepared implant site. The BTK dental implants can be positioned manually with the Reversible Torque Wrench or they can be inserted using the contra-angle handpiece. A range of 15 - 25 rpm is recommended for implant insertion and not to exceed the maximum torque indicated by BTK.

IMPLANT INSERTION



If a (partial / complete) tapping has been performed before implant insertion, the implant should be placed carefully on the implant site by performing a counterclockwise half rotation to engage the implant with the prepared thread.

After inserting into the thread, the implant can be guided in its final position clockwise, making sure that the implant is inserted at the desired depth and that the connection is intact.

WARNING The "implant driver" insertion drivers for the Morse-taper hexagon implant connection must be removed with delicate off-axis movements before removing them.



Insert the implant slowly in the previously prepared site.

A range of 15-25 rpm is recommended. During insertion, do not exceed the maximum torque values indicated below:

- **implant ≤ Ø 3,7 mm:**
Insertion torque max. 35 - 45 Ncm
- **implant > Ø 3,7 mm:**
Insertion torque max. 45 - 65 Ncm




In the cap of the internal vial there is, for each implant family, the corresponding cover screw (locking screw), sterile and ready for use.

Use sterile saline solution to carefully clean the implant connection from any organic residues. Therefore, make sure that it is clean and dry, before placing the cover screw (locking screw) or any prosthetic components that have been decided to be connect to the implant.

The cover screw is the chosen solution for the closed healing mode. To remove it more easily at the end of the healing period, a small amount of sterile vaseline or sterile chlorhexidine gel can be applied to the thread of the cover screw or healing cap before tightening it manually (5-8 Ncm) on the BTK implant, using a driver with a hex connection.



It is advisable to perform a postoperative x-ray check.







SIMPLICITY REDEFINED ONE KIT

The surgical tray is used for the secure storage and sterilization of surgical and auxiliary instruments.

The surgical tray is made of a highly shock-proof thermoplastic, which is well established in medical applications and the material is suitable for frequent sterilization in the autoclave.

General guidelines for the cleaning and sterilization are given in the corresponding GENERAL SURGICAL GUIDELINES (Ref. 06200117).

INSTRUMENTS KIT CODE 624NA001



1. GUIDE SHAFT
502MA002 | Ø 2.5 mm
2. ANGLED WRENCH 30°
502MA003 | HEX 3.10
3. DEPTH GAUGE
540MA011 | Ø 1.8 mm L108 mm 30°
4. SURGICAL GUIDE BT4
502MA006 | (PIN Ø2.5mm)
5. TORQUE WRENCH JD, REVERSIBLE
501JD003 | 90 Ncm



6

BONE PROFILER HS (BT4)

435HS430 | Ø 4.3 mm L25 mm

7

SCREWDRIVER JD (BT4)

530JD021 | HEX 2.0 L10 mm

530JD014 | HEX1.20 L15mm Slim Shank

8

PARALLELISM PIN (BT4)

540MA007 | M1.4 L26mm

9

BONE PROFILER GUIDE

435EN001 | EN
435ER001 | ER
435IR001 | IR
435KR001 | KR
435KW001 | KW

1

SHARP LANCE DRILL



401HR202 | Ø 2 mm L33 mm

2

TWIST DRILL



426HR200 | Ø 2 mm L36.5 mm

3

DRILLS



426HR250 | Ø 2,5 mm L36.5 mm



426HR310 | Ø 2,75 - 3,10 mm L36.5 mm



426HR345 | Ø 3,45 - 3,05 mm L36.5 mm



426HR385 | Ø 3,85 - 3,4 mm L36.5 mm



426HR420 | Ø 4.2 - 3.7 mm L36.8 mm



426HR455 | Ø 4.55 - 4 mm L36.8 mm



426HR540 | Ø 5.4 - 4.7 mm L36.8 mm



426HR570 | Ø 5.7 - 4.95 mm L36.8 mm

4

COUNTERSINK ISY KONE INT-EXT



434HS330 | L31 mm
Ø 3.3 mm (INT)
Ø 3.7 mm (INT/EXT)



434HS480 | L31 mm
Ø 4.8 mm (INT)

5

PARALLELISM PIN/ DEPTH GAUGE



540MA014 | Ø 2-2,5 mm
L36 mm 2 pcs



540MA019 | Ø 2-2,5 mm
L26 mm 2 pcs

6

TAPS ISY KONE



467HR290 | Ø 2,9 mm L36 mm



467HR330 | Ø 3,3 mm L36 mm



467HR370 | Ø 3,7 mm L36 mm



467HR410 | Ø 4,1 mm L36 mm



467HR480 | Ø 4,8 mm L36 mm



467HR600 | Ø 6 mm L36 mm


TAPS BT SAFE




466HR330 | Ø 3.3 mm L31 mm




466HR370 | Ø 3.7 mm L31 mm



466HR410 | Ø 4.1 mm L31 mm



466HR480 | Ø 4.8 mm L31 mm



466HR600 | Ø 6 mm L31 mm

7

IMPLANT DRIVER



530HL003 | KN (REGULAR) L33 mm



530HS013 | KR (SHORT) L23 mm



530HL001 | KR (REGULAR) L33 mm



530HS014 | KW (SHORT) L23 mm



530HL002 | KW (REGULAR) L=33 mm



530HS022 | IR-IW (SHORT) L25 mm



530HL005 | IR-IW (REGULAR) L35 mm



530HS023 | DR (SHORT) L25mm



530HL006 | DR (REGULAR) L35 mm

8

DRILLS EXTENSION



520HS003 | L28 mm

9

RETENTIVE WRENCH



530HS017 | HEX 3.10

10

SCREWDRIVER JD



530JD012 | HEX 0.90 L15 mm



530JD005 | HEX 1.20 L15 mm

11

MANUAL WRENCH - JD



530JD033 | ISO/HEX3.10-JD L10mm

12

DRILL STOP



518NA505 | H5 mm Snap Fit



518NA506 | H6 mm Snap Fit




518NA508 | H8 mm Snap Fit




518NA510 | H10 mm Snap Fit



518NA512 | H12 mm Snap Fit



518NA514 | H14 mm Snap Fit



518NA516 | H16 mm Snap Fit



518NA705 | H5 mm Snap Fit



518NA706 | H6 mm Snap Fit




518NA708 | H8 mm Snap Fit




518NA710 | H10 mm Snap Fit




518NA712 | H12 mm Snap Fit



518NA714 | H14 mm Snap Fit



518NA716 | H16 mm Snap Fit



518NA718 | H18 mm Snap Fit



518NA720 | H20 mm Snap Fit

CHARACTERISTICS OF SURGICAL DRILLS

- All drills and screw taps are made of stainless steel.
- All drills and screw taps are supplied in non-sterile single packs or in kit not sterile.
Please refer to the recommendations on cleansing and sterilization indicated by BTK.
- Drills and screw taps must be replaced after a maximum of 20 uses.
The effectiveness decreases after 5/6 applications already.
- All drills and screw taps have depth markings made with laser technique.
- The length relative to the corresponding black strip, realized with laser technique, it is always the lower or upper end of the strip.
- The black strips correspond to the length of the selected implant.
However, to increase security, the drill stops can be used during site preparation.
- All drills report their diameter and the relevant reference code on the stem.
- All final drills allow you to apply suitable drill stops.
- In case the length of the drills is insufficient, there is the possibility to connect them to the "Drill Extension" tool.

SURGICAL STANDARDS

For successful osseointegration, a precise, not traumatic surgical technique is required, which safeguards the soft tissues and accurately prepares the implant site without overheating the bone.


Before starting the surgical procedure and during the same procedure the following points must be taken into account:

- Check that all the necessary tools are available and fully functional.
It is recommended to always keep an adequate supply of sterile implants and instruments.
- Do not use cutting tools more than 20 times.
Make sure that the drills are sharp before each use.
The effectiveness of a drill already decreases after 5/6 applications.
- Drilling must be carried out with sharp drills, intermittently at 500 - 600 rpm, always with abundant external irrigation with pre-cooled sterile saline solution and avoiding excessive pressures.
- Do not exceed the speeds indicated by BTK for drills.
- Use the drills with diameters in ascending order.
- The drills can be placed in distilled / deionized water but should not be placed in saline or Ringer's solution during surgery if they are used for more than one preparation.

NOTE

- For implants with a length of 18 mm, the corresponding depth markings on drills are not provided such as the suitable drill stops are not available. It is advisable to prepare the implant site taking into account that the length of the drill, from the tip and up to the retention collar, is equal to 18.8 mm. It is responsibility of the clinician to evaluate based on the clinical case, morphology and bone quality, as well as the inclination of the implant, how to prepare the implant site.
- For ISY KONE implants with a diameter of Ø 6 mm, the corresponding screw tap (ref 467HR600 Screw Tap HR Ø6.0mm L36mm Isy Kone) is not included in the surgical kit (ref. 624NA001). This tool is OPTIONAL and must be purchased separately if necessary.

TOOL	SPEED (RPM)
LANCE DRILL	800
TWIST DRILL (PILOT) Ø 2mm	800
DRILL Ø<3.5mm	600
DRILL 3.5mm≤Ø≤4.5mm	500
DRILL Ø>4.5mm	400
COUNTERSINK	300/400
TAPS	<15




REVERSIBLE TORQUE WRENCH

The Reversible Torque Wrench is a dismantable, multiple-use instrument that provides means of tightening implants, abutments and screws. The lever arm integrated in the Reversible Torque Wrench is pushed away from the main body to the desired torque value. A torque value indicator is mounted at 90° in relation to the lever arm and indicates different value marks.

NOTE

Before the first and each following use, the Reversible Torque Wrench should be dismantled, cleaned, disinfected and sterilized in accordance with the instructions for use.



To dismantle the torque wrench for cleaning procedure, unscrew the wheel and then remove the inner bar where the spring is assembled.

APPLYING THE CORRECT TORQUE VALUE

In order to achieve the desired torque value, apply the force only to the lever-arm to the desired value mark. **The following marks are indicated: 15, 25, 35, 50, 70 and 90 Ncm.** Make sure that the arrow of the inversion device is matching to the lever-arm direction.

HOW TO CHANGE DIRECTIONS

With this type of Reversible Torque Wrench, one is able to change directions by simply pulling (1) and turning (2) the inversion device 180° in the desired direction.

This is done without removing the Reversible Torque Wrench from the attached driver in order to avoid additional manipulations and to save time.

The grey arrow on the inversion device always indicates in which direction the force is applied (3). This design was chosen to avoid additional manipulation, reduce potential sources of error while helping to save time.



DEVICE	MATERIAL	TIGHTENING TORQUE
Cover screw	Titanium GR5	from 5 to 8 Ncm (hand tight)
Healing abutment	Titanium GR5	from 5 to 8 Ncm (hand tight)
Impression Post screw, tightening to implant or implant replica	Titanium GR5	from 5 to 8 Ncm (hand tight)
Retentive screw, tightening Scan Abutment	Titanium GR5	from 5 to 8 Ncm (hand tight)
Retentive screw, temporary tightening (abutment to implant)	Titanium GR5	from 15 to 20 Ncm
Retentive screw, final tightening (abutment to implant)	Titanium GR5	from 25 to 30 Ncm
Straight abutment M.U.A.	Titanium GR5	from 25 to 30 Ncm
Retentive screw, tightening angled abutment M.U.A.	Titanium GR5	from 20 to 25 Ncm
Retentive screw, prosthesis to abutment M.U.A. - suprastructures	Titanium GR5	10 Ncm
Connector abutment to implant	Titanium GR5	from 20 to 25 Ncm

SURGICAL PROTOCOLS

DRILLS										TAPS (D1-D2)*												COUNTERSINKS (D1-D2)*					
IMPLANT		Lance drill 401HR202 ***	Twist drill (Pilot) Ø 2,00 426HR200 ***	Twist drill Ø 2,5 426HR250	Twist drill Ø 3,1 426HR310	Twist drill Ø 3,45 426HR345	Twist drill Ø 3,85 426HR385	Twist drill Ø 4,2 426HR420	Twist drill Ø 4,55 426HR455	Twist drill Ø 5,4 426HR540	Twist drill Ø 5,7 426HR570	Screw tap Ø2,9 467HR290		Screw tap Ø3,3 467HR330		Screw tap Ø3,7 467HR370		Screw tap Ø4,1 467HR410		Screw tap Ø4,8 467HR480		Screw tap Ø6 467HR600		Countersink Ø3,3/ Ø3,7 434HS330		Countersink Ø4,8 434HS480	
												D1	D2	D1	D2	D1	D2	D1	D2	D1	D2	D1	D2	D1	D2	D1	D2
2,9	8mm											-	-												(IR only)		
	10mm	●	●	○								-	-														
	12mm											1°**	1°**														
	14mm											1°**	1°**														
3,3	8mm													1°	-										(IR only)		
	10mm	●	●	●	○									1°	-												
	12mm													1°	1°												
	14mm													1°	1°												
3,7	8mm															1°	-								(IR and ER)		
	10mm															1°	-										
	12mm	●	●	●	●	○										1°	1°										
	14mm															1°	1°										
	16mm															2°	1°										
	18mm															2°	1°										
4,1	8mm																	1°	1°								
	10mm																	2°	1°								
	12mm	●	●	●	●	●	○											3°	2°								
	14mm																	3°	2°								
	16mm																	3°	3°								
	18mm																	3°	3°								
4,8	8mm																			1°	1°				(only IW)		
	10mm	●	●	●	●	●	●	○												2°	1°						
	12mm																			3°	2°						
	14mm																			3°	2°						
6	8mm	●	●	●	●	●	●	●	●	●	○											1°	1°				
	10mm																					1°	1°				

● Always

○ Only in presence of D1 bone

○ Only in presence of D2 and D3 bone

** No Screw taps and Countersinks are necessary for D3 and D4

*** Screw taps should be used when the standard protocol is not sufficient to insert completely the implant without oltrepassing the maximum insertion torque recommended (max 45 Ncm) and it is a choice that must be evaluated by the clinician after a careful analysis


*** This tool is included in the surgical kit ref, 624NA001

NOTE: ISY KONE 2.9 is not suggested for D4 bone

N.B.

The above procedures should be considered indicative; it is responsibility of the clinician to evaluate potential variations of the procedure on the basis of individual case and bone density.





























Due to the presence of cutting edges on drills, it is strongly recommended the use of the stop in order to avoid excessive depth of drilling that could compromise vital structures.


































ADDITIONAL INSTRUMENTS












Optional instrumentation dedicated to the correct management of the surgical procedure.

PICTURE	REF	PRODUCT NAME	SPECIFICATION
INITIAL PREPARATION			
	401HS200	Round Drill HS	Ø2mm L30mm
	401HS201	Lance Drill HS	Ø2mm L30mm
	401HR200	Round Drill HR	Ø2mm L35mm
	401HR201	Lance Drill HR	Ø2mm L35mm
DRILLS, LENGHT 32.5-32.8 mm (SHORT)			
	426HS200	Twist Drill HS	Ø2mm L32.5mm
	426HS250	Twist Drill HS	Ø2.5mm L32.5mm
	426HS270	Twist Drill HS	Ø2.7mm L32.5mm
	426HS310	Twist Step Drill HS	Ø3.1-2.75mm L32.5mm
	426HS345	Twist Step Drill HS	Ø3.45-3.05mm L32.5mm
	426HS385	Twist Step Drill HS	Ø3.85-3.4mm L32.8mm
	426HS420	Twist Step Drill HS	Ø4.2-3.7mm L32.8mm
	426HS455	Twist Step Drill HS	Ø4.55-4mm L32.8mm
	426HS540	Twist Step Drill HS	Ø5.4-4.7mm L32.8mm
	426HS570	Twist Step Drill HS	Ø5.7-4.95mm L32.8mm
DRILL STOPS (SNAP-FIT) Ø5 FOR SHORT DRILLS ≤ Ø3.45 mm			
	690NA258	Stop Kit	Ø5 S5-12mm BT Safe Isy Kone BT Nano
	521NA505	Drill Stop	H5mm Snap Fit Short
	521NA506	Drill Stop	H6mm Snap Fit Short
	521NA507	Drill Stop	H7mm Snap Fit Short
	521NA508	Drill Stop	H8mm Snap Fit Short
	521NA510	Drill Stop	H10mm Snap Fit Short
	521NA512	Drill Stop	H12mm Snap Fit Short
DRILL STOPS (SNAP-FIT) Ø7 FOR SHORT DRILLS ≥ Ø3.85 mm			
	690NA259	Stop Kit	Ø7 S5-12mm BT Safe Isy Kone BT Nano
	521NA705	Drill Stop	H5mm Snap Fit Short
	521NA706	Drill Stop	H6mm Snap Fit Short
	521NA707	Drill Stop	H7mm Snap Fit Short
	521NA708	Drill Stop	H8mm Snap Fit Short
	521NA710	Drill Stop	H10mm Snap Fit Short
	521NA712	Drill Stop	H12mm Snap Fit Short














PICTURE	REF	PRODUCT NAME	SPECIFICATION	
DRILL STOPS (SNAP-FIT) Ø5 FOR REGULAR DRILLS ≤ Ø3.45 mm				
	690NA256	Kit Stop	Ø5 R5-16mm BT Safe Isy Kone BT Nano	
DRILL STOPS (SNAP-FIT) Ø7 FOR REGULAR DRILLS ≥ Ø3.85 mm				
	690NA257	Kit Stop	Ø7 R5-16mm BT Safe Isy Kone BT Nano	
TAP DRIVERS FOR MANUAL USE COMPATIBLE WITH REVERSIBLE TORQUE WRENCH (JD)				
	530JD031	Adapter Connection	ISO/HEX3.10-JD L35mm	
	530JD032	Adapter Connection	ISO/HEX3.10-JD L7.5mm	
	530JD034	Adapter Connection	ISO/HEX3.10-JD L15mm	
MOUNTING DEVICES (FOR ISY KONE EN/ER)				
EN		690EN001	Mounting Device EN	L8mm
		690EN002	Mounting Device EN	L20mm
ER		690ER001	Mounting Device ER	L8mm
		690ER002	Mounting Device ER	L20mm
IMPLANT DRIVERS, FOR MANUAL USE COMPATIBLE WITH TORQUE WRENCH (JD)				
KR		530JD023	Implant Driver JD Conn. KR	L15mm
		530JD024	Implant Driver JD Conn. KR	L23mm
KW		530JD027	Implant Driver JD Conn. KW	L15mm
		530JD028	Implant Driver JD Conn. KW	L23mm
IR-IW		530JD018	Implant Driver JD Conn. IR-IW	L19mm
		530JD019	Implant Driver JD Conn. IR-IW	L22mm
HANDPIECE DRIVER				
	530HS002	Handpiece Driver	HEX0.90 L25mm	
	530HS003	Handpiece Driver	HEX0.90 L30mm	
	530HS004	Handpiece Driver	HEX1.20 L25mm	
	530HS005	Handpiece Driver	HEX1.20 L30mm	
HEX DRIVERS FOR MANUAL USE COMPATIBLE WITH REVERSIBLE TORQUE RATCHET (JD)				
	530JD003	Screwdriver JD	HEX1.20 L5mm	
	530JD004	Screwdriver JD	HEX1.20 L10mm	
	530JD006	Screwdriver JD	HEX1.20 L20mm	
	530JD007	Screwdriver JD	HEX1.20 L30mm	
	530JD011	Screwdriver JD	HEX0.90 L10mm	
INSTRUMENTS FOR BT4 METHOD				
	530HS012	Handpiece Driver	HEX1.20 L30mm Reduced	
	530JD015	Screwdriver JD	HEX2.0 L5mm	
	530JD038	Screwdriver JD	HEX2.0 L20mm	
	530JD036	Screwdriver JD	HEX1.50 L15mm	
	530JD037	Screwdriver JD	HEX1.50 L30mm	

MORSE-TAPER (MTH) KN/KR





HEALING & SOFT TISSUE CONDITIONING

KN	KR	PICTURE	REF	PRODUCT NAME	SPECIFICATION
HEALING ABUTMENTS Ø 1.8 mm					
			201KR5A3	Healing Abutments KR	H5mm Ø1.8mm
			201KR7A0	Healing Abutments KR	H7mm Ø1.8mm
HEALING ABUTMENTS Ø 2.5 mm					
○			201KN2A2	Healing Abutments KN	H2mm Ø2.5mm
			201KN3A2	Healing Abutments KN	H3.5mm Ø2.5mm
			201KN5A2	Healing Abutments KN	H5mm Ø2.5mm
	○		201KR2A5	Healing Abutments KR	H2mm Ø2.5mm
			201KR3A4	Healing Abutments KR	H3.5mm Ø2.5mm
			201KR5A1	Healing Abutments KR	H5mm Ø2.5mm
HEALING ABUTMENTS Ø 3.5 mm					
○			201KN2A0	Healing Abutments KN	H2mm Ø3.5mm
			201KN3A0	Healing Abutments KN	H3.5mm Ø3.5mm
	○		201KR1A0	Healing Abutments KR	H1mm Ø3.5mm
			201KR2A1	Healing Abutments KR	H2mm Ø3.5mm
			201KR3A0	Healing Abutments KR	H3.5mm Ø3.5mm
			201KR5A4	Healing Abutments KR	H5mm Ø3.5mm
			201KR7A2	Healing Abutments KR	H7mm Ø3.5mm
HEALING ABUTMENTS Ø 4 mm					
	○		201KR6A0	Healing Abutments KR	H6mm Ø4mm
HEALING ABUTMENTS Ø 4.5 mm					
○			201KN2A1	Healing Abutments KN	H2mm Ø4.5mm
			201KN3A1	Healing Abutments KN	H3.5mm Ø4.5mm
	○		201KR1A1	Healing Abutments KR	H1mm Ø4.5mm
			201KR2A2	Healing Abutments KR	H2mm Ø4.5mm
			201KR3A1	Healing Abutments KR	H3.5mm Ø4.5mm
			201KR5A2	Healing Abutments KR	H5mm Ø4.5mm
			201KR7A1	Healing Abutments KR	H7mm Ø4.5mm
HEALING ABUTMENTS Ø 5.5 mm					
	○		201KR2A3	Healing Abutments KR	H2mm Ø5.5mm
			201KR3A2	Healing Abutments KR	H3.5mm Ø5.5mm
			201KR5A5	Healing Abutments KR	H5mm Ø5.5mm
HEALING ABUTMENTS Ø 6.5 mm					
	○		201KR2A4	Healing Abutments KR	H2mm Ø6.5mm
			201KR3A3	Healing Abutments KR	H3.5mm Ø6.5mm
			201KR5A6	Healing Abutments KR	H5mm Ø6.5mm
MILLABLE HEALING ABUTMENTS					
	○		204KR0A0	Millable Healing Abutments KR	H4mm Ø4.3mm
			204KR0A1	Millable Healing Abutments KR	H4mm Ø5.5mm
			204KR0A2	Millable Healing Abutments KR	H4mm Ø6mm
			204KR0A3	Millable Healing Abutments KR	H4mm Ø7.2mm

















IMPRESSION TAKING

KN	KR	PICTURE	REF	PRODUCT NAME	SPECIFICATION
IMPRESSION POST					
	○		325KR0A0	Impression Post Pro KR	Plastic cap
			690NA073	Impression Post Screw	M1.6 HEX1.20 H7.9mm
			690NA091.10	Caps Kit Pro	Ø5.1mm Kit 10pz
IMPRESSION POST PROPICK-UP					
○			323KN0A0	Transfer Propick-Up KN	HUseful21.5mm Long Screw
○			690NA299	Impression Post Pick-Up Screw	M1.4 HEX1.20 H26mm
	○		323KR0A0	Impression Post Propick-Up KR	HUseful21.5mm Long Screw
			323KR0A2	Impression Post Propick-Up KR	HUseful16.5mm Short Screw
			323KR0R0	Impression Post Propick-Up KR	HUseful21.5mm Rotating Long Screw
			323KR0R1	Impression Post Propick-Up KR	HUseful16.5mm Rotating Short Screw
	○		690NA071	Impression Post Pick-Up Screw	M1.6 HEX1.20 H21.4mm
	○		690NA072	Impression Post Pick-Up Screw	M1.6 HEX1.20 H26.4mm
IMPLANT REPLICA					
○			301KN0A0	Implant Replica KN	
	○		301KR0A0	Implant Replica KR	

INTERIM RESTORATIONS

KN	KR	PICTURE	REF	PRODUCT NAME	SPECIFICATION
TEMPORARY ABUTMENTS					
○			210KN1A0	Temporary Abutment KN	
	○		210KR1A0	Temporary Abutment KR	
			210KR1R0	Temporary Abutment KR	Rotating
			215KR0A0	Temporary Abutment KR	Peek

CEMENT-RETAINED PROSTHESIS

KN	KR	PICTURE	REF	PRODUCT NAME	SPECIFICATION
	○		214KR4A0	Transfer Abutment KR	H4mm Ø4.5mm
	○		690NA091.10	Caps Kit Pro	Ø5.1mm Kit 10pcs
STRAIGHT AND ESTHETIC ABUTMENTS Ø 3.5 mm					
	○		220KR1A2	Straight Abutment KR	H1mm Ø3.5mm
			220KR2A0	Straight Abutment KR	H2mm Ø3.5mm
			220KR3A0	Straight Abutment KR	H3.5mm Ø3.5mm
○			219KN2A0	Esthetic Abutment KN	H2mm Ø3.5mm
			219KN3A0	Esthetic Abutment KN	H3.5mm Ø3.5mm
STRAIGHT AND ESTHETIC ABUTMENTS Ø 4.5 mm					
	○		220KR1A3	Straight Abutment KR	H1mm Ø4.5mm
			220KR2A1	Straight Abutment KR	H2mm Ø4.5mm
			220KR3A1	Straight Abutment KR	H3.5mm Ø4.5mm
○			219KN2A1	Esthetic Abutment KN	H2mm Ø4.5mm
			219KN3A1	Esthetic Abutment KN	H3.5mm Ø4.5mm
	○		219KR2A3	Esthetic Abutment KR	H2mm Ø4.5mm
			219KR3A3	Esthetic Abutment KR	H3.5mm Ø4.5mm
STRAIGHT AND ESTHETIC ABUTMENTS Ø 5.5 mm					
	○		220KR1A1	Straight Abutment KR	H1.5mm Ø5.5mm
			220KR2A2	Straight Abutment KR	H2mm Ø5.5mm
			220KR3A2	Straight Abutment KR	H3.5mm Ø5.5mm
	○		219KR2A4	Esthetic Abutment KR	H2mm Ø5.5mm
			219KR3A1	Esthetic Abutment KR	H3mm Ø5.5mm
			219KR3A4	Esthetic Abutment KR	H3.5mm Ø5.5mm
ESTHETIC ABUTMENTS Ø 6.5 mm					
	○		219KR2A2	Esthetic Abutment KR	H2mm Ø6.5mm
			219KR3A2	Esthetic Abutment KR	H3.5mm Ø6.5mm
ANGLED ABUTMENT 10°					
○			220KN2C0	Angled Abutment KN	10° H2mm Ø3.5mm
			220KN2C1	Angled Abutment KN	10° H2mm Ø4.5mm
			220KN3C0	Angled Abutment KN	10° H3.5mm Ø3.5mm
			220KN3C1	Angled Abutment KN	10° H3.5mm Ø4.5mm
	○		220KR1C0	Angled Abutment KR	10° H1mm Ø3.5mm
			220KR1C1	Angled Abutment KR	10° H1mm Ø4.5mm
			220KR2C0	Angled Abutment KR	10° H2mm Ø3.5mm
			220KR2C1	Angled Abutment KR	10° H2mm Ø4.5mm
			220KR2C2	Angled Abutment KR	10° H2mm Ø5.5mm
			220KR3C0	Angled Abutment KR	10° H3.5mm Ø3.5mm
			220KR3C1	Angled Abutment KR	10° H3.5mm Ø4.5mm
			220KR3C2	Angled Abutment KR	10° H3.5mm Ø5.5mm
ANGLED ABUTMENT 20°					
	○		220KR1E0	Angled Abutment KR	20° H1mm Ø3.5mm
			220KR1E1	Angled Abutment KR	20° H1mm Ø4.5mm
			220KR2E2	Angled Abutment KR	20° H2mm Ø3.5mm
			220KR2E3	Angled Abutment KR	20° H2mm Ø4.5mm
			220KR2E4	Angled Abutment KR	20° H2mm Ø5.5mm
			220KR3E0	Angled Abutment KR	20° H3.5mm Ø3.5mm
			220KR3E1	Angled Abutment KR	20° H3.5mm Ø4.5mm
			220KR3E2	Angled Abutment KR	20° H3.5mm Ø5.5mm
			220KR4E1	Angled Abutment KR	20° H4mm Ø5.5mm
ANGLED ABUTMENT 30°					
	○		220KR1G0	Angled Abutment KR	30° H1mm Ø3.5mm
			220KR1G1	Angled Abutment KR	30° H1mm Ø4.5mm
			220KR2G0	Angled Abutment KR	30° H2mm Ø3.5mm
			220KR2G1	Angled Abutment KR	30° H2mm Ø4.5mm
			220KR2G2	Angled Abutment KR	30° H2mm Ø5.5mm
			220KR3G0	Angled Abutment KR	30° H3.5mm Ø3.5mm
			220KR3G1	Angled Abutment KR	30° H3.5mm Ø4.5mm
			220KR3G2	Angled Abutment KR	30° H3.5mm Ø5.5mm
RETENTIVE SCREWS FOR ABUTMENTS					
○			690NA297	Retentive Screw KN	M1.4 HEX1.20 H8mm
	○		690NA070	Retentive Screw KR	M1.6 HEX1.20 H8.3mm

SCREW-RETAINED / CEMENT-RETAINED PROSTHESIS













KN	KR	PICTURE	REF	PRODUCT NAME	SPECIFICATION
BT LINK					
○			246KN1A0	BT Link KN	H1mm Ø4.1mm
			246KN1A2	BT Link KN	H1mm Ø3.6mm
○			247KN1A0	Base BT Link KN	H1mm Ø4.1mm no Cap.
			247KN1A2	Base BT Link KN	H1mm Ø3.6mm no Cap.
○			690NA298	Retentive Screw KN	M1.4 HEX1.20 H8mm TP
	○		246KR1A0	BT Link KR	H1mm Ø4.1mm
			246KR1A1	BT Link KR	H1mm Ø4.1mm Rotating
	○		247KR1A0	Base BT Link KR	H1mm Ø4.1mm no Cap.
			247KR1A1	Base BT Link KR	H1mm Ø4.1mm no Cap. Rot.
	○		244KR0A0	SIR Link KR	H0.8mm Ø4.1mm
	○		690NA083	Retentive Screw KR	M1.6 HEX1.20 H8.3mm TP
○	○		205NA001.05	Castable Plastic Abut. BT Link	H1mm Ø4.7mm Kit 5pcs
BT GRIP					
	○		530JD036	Screwdriver JD BT GRIP	HEX1.50 L15 mm (Short)
	○		530JD037	Screwdriver JD BT GRIP	HEX1.50 L30 mm (Long)
	○		248ER1A0	X3 Link KR	H1mm Ø4.1mm Multi High
	○		690NA239	Retentive Screw BT GRIP	M1.6 HEX1.50
CAST-ON TECHNIQUE					
○			240KN1A0	CoCr Abutment KN	H1.5mm
○			690NA297	Retentive Screw KN	M1.4 HEX1.20 H8mm
	○		245KR0A0	Gold Abutment KR	H1mm
			240KR1A0	CoCr Abutment KR	H1.5mm
			240KR1R0	CoCr Abutment KR	H1.5mm Rotating
	○		690NA070	Retentive Screw KR	M1.6 HEX1.20 H8.3mm
SCAN ABUTMENT					
○			351KN1A0	Scan Abutment Extra-oral KN	
○			352KN1A0	Scan Abutment Intra-oral KN	
○			690NA298	Retentive Screw KN	M1.4 HEX1.20 H8mm TP
	○		351KR1A0	Scan Abutment Extra-oral KR	
	○		352KR1A0	Scan Abutment Intra-oral KR	
	○		690NA083	Retentive Screw KR	M1.6 HEX1.20 H8.3mm TP

SCREW-RETAINED PROSTHESIS




















KR	PICTURE	REF	PRODUCT NAME	SPECIFICATION
BT4 STRAIGHT ABUTMENTS				
		265KR1R0	BT4 Straight Abutment KR	Rotating H1mm
		265KR2R0	BT4 Straight Abutment KR	Rotating H2mm
		265KR3R0	BT4 Straight Abutment KR	Rotating H3mm
		268KR1R0	BT4 Slim Straight Abutment KR	Rotating H1mm
		268KR2R0	BT4 Slim Straight Abutment KR	Rotating H2mm
BT4 ANGLED ABUTMENTS				
		266KR2L0	BT4 Angled Abutment KR	17° H2mm Ø4.8mm
		266KR4L0	BT4 Angled Abutment KR	17° H4mm Ø4.8mm
		266KR3L0	BT4 Angled Abutment KR	17° H3mm Ø 4.8mm
		266KR3G0	BT4 Angled Abutment KR	30° H3mm Ø 4.8mm
		690NA075	Retentive Screw BT4	M1.6 Angled Abutment KR
COVERING CAPS BT4				
		330NA0A0.04	Covering Caps BT4	H5 Kit 4pcs
		330BU0A0.04	Covering Caps BT4 Slim	H5 Kit 4pcs
		690NA024	Retentive Screw	M1.4 HEX1.20 10N
BT4 INTERIM RESTORATION				
		267NA0A0	BT4 Titanium Abutment	
		269BU0A0	BT4 Slim Titanium Abutment	
		207NA0A0	Castable Plastic Abutment BT4	
		207NA0A1	Castable Plastic Abutment BT4	No screw
		207BU1R0	Castable Plastic Abutment BT4 Slim	
		690NA024	Retentive Screw	M1.4 HEX1.20 10N
		311NA0A0	Impression Post Pick-up BT4	with long Screw
		311BU0A0	Impression Post Pick-up BT4 Slim	with long Screw
		690NA031	Impression Post Pick-Up Screw	M1.4 HEX1.20 H17mm
		303NA0A0	Abutment replica BT4	
		303BU0A0	Abutment replica BT4 Slim	
BT4 SCAN ABUTMENT				
		351BT1A1	Scan Abutment Extra-oral BT	Rotating
		352BT1A1	Scan Abutment Intra-oral BT	Rotating
		351BU1A1	Scan Abutment Extra-oral BU	Rotating (BT4 Slim)
		352BU1A1	Scan Abutment Intra-oral BU	Rotating (BT4 Slim)
		690NA024	Retentive Screw	M1.4 HEX1.20 10N
BT4 BT LINK				
		246BT1A1	BT Link BT	H1mm Ø4.8mm Rotating
		247BT1A1	Base BT Link BT	H1mm Ø4.8mm no Cap. Rot.
		205NA003.05	Castable Plastic Abut. BT Link	H1mm Ø5.4mm Kit 5pcs
		246BU1A1	BT Link BU	H1mm Ø4.1mm Rotating
		247BU1A1	Base BT Link BU	H1mm Ø4.1mm no Cap. Rot.
		205NA001.05	Castable Plastic Abut. BT Link	H1mm Ø4.7mm Kit 5pcs
		690NA024	Retentive Screw	M1.4 HEX1.20 10N
CAST-ON TECHNIQUE				
		240BT1R0	Base CoCr BT	H1.5mm Rotating
		240BU1R0	Base CoCr BU	H1.5mm Rotating
		690NA024	Retentive Screw	M1.4 HEX1.20 10N

INTERNAL HEXAGON (INT) IR/IW






HEALING & SOFT TISSUE CONDITIONING

IR	IW	PICTURE	REF	PRODUCT NAME	SPECIFICATION
HEALING ABUTMENTS Ø 4.5 mm					
			201IR2A0	Healing Abutment IR	H2mm Ø4.5mm
			201IR3A0	Healing Abutment IR	H3mm Ø4.5mm
			201IR4A0	Healing Abutment IR	H4mm Ø4.5mm
			201IR5A0	Healing Abutment IR	H5mm Ø4.5mm
HEALING ABUTMENTS Ø 5.5 mm					
			201IR4A1	Healing Abutment IR	H4mm Ø5.5mm
			201IW5A0	Healing Abutment IW	H5mm Ø5.5mm
HEALING ABUTMENTS Ø 6 mm					
			201IW2A0	Healing Abutment IW	H2mm Ø6mm
			201IW3A0	Healing Abutment IW	H3mm Ø6mm
			201IW4A0	Healing Abutment IW	H4mm Ø6mm































IMPRESSION TAKING

IR	IW	PICTURE	REF	PRODUCT NAME	SPECIFICATION
IMPRESSION POST					
●			320IR0A1	Impression Post IR	Aluminium Cap
	●		320IW0A1	Impression Post IW	Aluminium Cap
●	●		690NA029	Impression Post Caps	Aluminium Kit 3pcs
●			325IR0A0	Impression Post Pro IR	Plastic Cap
	●		325IW0A1	Impression Post Pro IW	Plastic Cap
●	●		690NA091.10	Caps Kit Pro	Ø5.1mm Kit 10pcs
●	●		690NA068	Impression Post Screw	M1.8 HEX1.20 H7.9mm
IMPRESSION POST PROPICK-UP					
●			323IR0A0	Impression Post Propick-Up IR	HUSeul21.5mm Long screw
			323IR0A1	Impression Post Propick-Up IR	HUSeul16.5mm Short screw
			323IR0R0	Impression Post Propick-Up IR	HUSeul21.5mm Rotating Long screw
			323IR0R1	Impression Post Propick-Up IR	HUSeul16.5mm Rotating Short screw
	●		323IW0A0	Impression Post Propick-Up IW	HUSeul21.5mm Long screw
			323IW0A1	Impression Post Propick-Up IW	HUSeul16.5mm Short screw
			323IW0R0	Impression Post Propick-Up IW	HUSeul21.5mm Rotating Long screw
			323IW0R1	Impression Post Propick-Up IW	HUSeul16.5mm Rotating Short screw
●	●		690NA222	Impression Post Pick-Up Screw	M1.8 HEX1.20 H26mm
			690NA190	Impression Post Pick-Up Screw	M1.8 HEX1.20 H21mm
IMPLANT REPLICA					
●			301IR0A1	Implant Replica IR	
	●		301IW0A1	Implant Replica IW	





























INTERIM RESTORATIONS

IR	IW	PICTURE	REF	PRODUCT NAME	SPECIFICATION
TEMPORARY ABUTMENTS					
●			210IR2A0	Temporary Abutment IR	
			210IR2A1	Temporary Abutment IR	Rotating
●			215IR2A0	Temporary Abutment IR	Peek
	●		210IW2A0	Temporary Abutment IW	
			210IW2R0	Temporary Abutment IW	Rotating
	●		215IW2A0	Temporary Abutment IW	Peek
●	●		690NA012	Retentive Screw	M1.8 HEX1.20






































CEMENT-RETAINED PROSTHESIS

IR	IW	PICTURE	REF	PRODUCT NAME	SPECIFICATION
STRAIGHT ABUTMENTS Ø 4.5 mm					
			220IR0A0	Straight Abutment IR	H0mm Ø4.5mm
			220IR2A1	Straight Abutment IR	H2mm Ø4.5mm
			220IR3A0	Straight Abutment IR	H3mm Ø4.5mm
			220IR4A0	Straight Abutment IR	H4mm Ø4.5mm
			220IR5A0	Straight Abutment IR	H5mm Ø4.5mm
STRAIGHT ABUTMENTS Ø 5.5 mm					
			220IW0A0	Straight Abutment IW	H0mm Ø5.5mm
			220IW2A0	Straight Abutment IW	H2mm Ø5.5mm
			220IW3A0	Straight Abutment IW	H3mm Ø5.5mm
			220IW5A0	Straight Abutment IW	H5mm Ø5.5mm
STRAIGHT ABUTMENTS Ø 6 mm					
			220IW4A0	Straight Abutment IW	H4mm Ø6mm
AESTHETIC ABUTMENTS					
			219IR0A0	Aesthetic Abutment IR	
			219IW0A0	Aesthetic Abutment IW	
ANGLED ABUTMENT 15°					
			220IR2D0	Angled Abutment IR	15° H2mm
			220IR4D0	Angled Abutment IR	15° H4mm
			220IW2D0	Angled Abutment IW	15° H2mm
ANGLED ABUTMENT 25°					
			220IR2F0	Angled Abutment IR	25° H2mm
			220IR4F0	Angled Abutment IR	25° H4mm
RETENTIVE SCREWS FOR ABUTMENT					
			690NA012	Retentive Screw	M1.8 HEX1.20
			VLE14TIT	Lingual Screw Esthetic Abut.	M1.4 HEX0.9

SCREW-RETAINED / CEMENT-RETAINED PROSTHESIS
















IR	IW	PICTURE	REF	PRODUCT NAME	SPECIFICATION
BT LINK					
●			246IR1A0	BT Link IR	H1mm Ø4.1mm
			246IR1A1	BT Link IR	H1mm Ø4.1mm Rotating
●			247IR1A0	Base BT Link IR	H1mm Ø4.1mm no Cap.
			247IR1A1	Base BT Link IR	H1mm Ø4.1mm no Cap. Rot.
●			205NA001.05	Castable Plastic Abut. BT Link	H1mm Ø4.7mm Kit 5pcs
●			244IR0A0	SIR Link IR	H0.8mm Ø4.1mm
	●		246IW1A0	BT Link IW	H1mm Ø5.6mm
			246IW1A1	BT Link IW	H1mm Ø5.6mm Rotating
	●		247IW1A0	Base BT Link IW	H1mm Ø5.6mm no Cap.
			247IW1A1	Base BT Link IW	H1mm Ø5.6mm no Cap. Rot.
	●		205NA005.05	Castable Plastic Abut. BT Link	H1mm Ø5.6mm Kit 5pcs
●	●		690NA077	Retentive Screw	M1.8 HEX1.20 H7.8mm TP
BT GRIP					
●			530JD036	Screwdriver JD BT GRIP	HEX1.50 L15 mm (Short)
			530JD037	Screwdriver JD BT GRIP	HEX1.50 L30 mm (Long)
●			248IR1A0	X3 Link IR	H1mm Ø4.1mm Multi High
●			690NA237	Retentive Screw BT GRIP	M1.8 HEX 1.50
CAST-ON TECHNIQUE					
●			245IR1A0	Gold Abutment IR	H1mm
			245IR1R0	Gold Abutment IR	H1mm Rotating
●			240IR1A0	CoCr Abutment IR	H1.5mm
			240IR1R0	CoCr Abutment IR	H1.5mm Rotating
	●		245IW1A0	Gold Abutment IW	H1mm
			245IW1R0	Gold Abutment IW	H1mm Rotating
	●		240IW1A0	CoCr Abutment IW	H1.5mm
			240IW1R0	CoCr Abutment IW	H1.5mm Rotating
●	●		690NA012	Retentive Screw	M1.8 HEX1.20
●	●		690NA018	Retentive Screw Gold	M1.8 HEX1.20
CASTABLE PLASTIC ABUTMENTS					
●			205IR1A0	Castable Plastic Abutment IR	
			205IR1A0.10	Castable Plastic Abutment IR	Kit 10pcs
●			205IR1A1	Castable Plastic Abutment IR	Rotating
			205IR1A1.10	Castable Plastic Abutment IR	Rotating Kit 10pcs
	●		205IW1A0	Castable Plastic Abutment IW	
			205IW1A0.10	Castable Plastic Abutment IW	Kit 10pcs
	●		205IW1R0	Castable Plastic Abutment IW	Rotating
			205IW1R0.10	Castable Plastic Abutment IW	Rotating Kit 10pcs
●	●		690NA012	Retentive Screw	M1.8 HEX1.20
●	●		690NA018	Retentive Screw Gold	M1.8 HEX1.20
SCAN ABUTMENT					
●			351IR1A0	Scan Abutment Extra-oral IR	
			352IR1A0	Scan Abutment Intra-oral IR	
	●		351IW1A0	Scan Abutment Extra-oral IW	
			352IW1A0	Scan Abutment Intra-oral IW	
●	●		690NA077	Retentive Screw	M1.8 HEX1.20 H7.8mm TP

SCREW-RETAINED PROSTHESIS

























IR	PICTURE	REF	PRODUCT NAME	SPECIFICATION
BT4 STRAIGHT ABUTMENT				
		265IR1R0	BT4 Straight Abutment IR	Rotating H1mm
		265IR2R0	BT4 Straight Abutment IR	Rotating H2mm
		265IR3R0	BT4 Straight Abutment IR	Rotating H3mm
ANGLED ABUTMENT BT4 17°				
		266IR2L0	BT4 Angled Abutment IR	17° H2mm Ø4.8mm
		266IR3L0	BT4 Angled Abutment IR	17° H3mm Ø4.8mm
ANGLED ABUTMENT BT4 30°				
		266IR3G0	BT4 Angled Abutment IR	30° H3mm Ø4.8mm
ANGLED ABUTMENT BT4 45°				
		266IR3M0	BT4 Angled Abutment IR	45° H3.5mm Ø5.7mm
ANGLED ABUTMENT BT4 55°				
		266IR3J0	BT4 Angled Abutment IR	55° H3.5mm Ø6.2mm
RETENTIVE SCREW BT4 FOR ANGLED ABUTMENT				
		690NA066	Retentive Screw BT4	M1.8 Angled abutment IR-EN
BT4 CAPS				
		330NA0A0.04	Covering Caps BT4	H5 Kit 4pcs
		690NA024	Retentive Screw	M1.4 HEX1.20 10N
BT4 INTERIM RESTORATION				
		267NA0A0	BT4 Titanium Abutment	
		207NA0A0	Castable Plastic Abutment BT4	
		207NA0A1	Castable Plastic Abutment BT4	No screw
		690NA024	Retentive Screw	M1.4 HEX1.20 10N
		311NA0A0	Impression Post Pick-up BT4	with Long Screw
		690NA031	Impression Post Pick-Up Screw	M1.4 HEX1.20 H17mm
		303NA0A0	Abutment Replica BT4	
BT4 SCAN ABUTMENT				
		351BT1A1	Scan Abutment Extra-oral BT	Rotating
		352BT1A1	Scan Abutment Intra-oral BT	Rotating
		690NA024	Retentive Screw	M1.4 HEX1.20 10N
BT4 BT LINK				
		246BT1A1	BT Link BT	H1mm Ø4.8mm Rotating
		247BT1A1	Base BT Link BT	H1mm Ø4.8mm no Cap. Rot
		205NA003.05	Castable Plastic Abut. BT Link	H1mm Ø5.4mm Kit 5pcs
		690NA024	Retentive Screw	M1.4 HEX1.20 10N
CAST-ON TECHNIQUE				
		240BT1R0	Base CoCr BT	H1.5mm Rotating
		690NA024	Retentive Screw	M1.4 HEX1.20 10N















EXTERNAL HEXAGON (EXT) EN/ER

HEALING & SOFT TISSUE CONDITIONING







EN	ER	PICTURE	REF	PRODUCT NAME	SPECIFICATION
HEALING ABUTMENTS Ø 4.5 mm					
			201EN2A0	Healing Abutment EN	H2mm Ø4.5mm
			201EN3A0	Healing Abutment EN	H3mm Ø4.5mm
			201EN4A0	Healing Abutment EN	H4mm Ø4.5mm
			201EN6A0	Healing Abutment EN	H6mm Ø4.5mm
HEALING ABUTMENTS Ø 5 mm					
			201ER2A0	Healing Abutment ER	H2mm Ø5mm
			201ER3A0	Healing Abutment ER	H3mm Ø5mm
			201ER4A0	Healing Abutment ER	H4mm Ø5mm
			201ER6A0	Healing Abutment ER	H6mm Ø5mm
			201ER8A0	Healing Abutment ER	H8mm Ø5mm
HEALING ABUTMENTS Ø 6 mm					
			201ER4A1	Healing Abutment ER	H4mm Ø6mm
			201ER6A1	Healing Abutment ER	H6mm Ø6mm
HEALING ABUTMENTS Ø 7 mm					
			201ER6A3	Healing Abutment ER	H6mm Ø7mm
HEALING ABUTMENTS Ø 7.5 mm					
			201ER4A2	Healing Abutment ER	H4mm Ø7.5mm

IMPRESSION TAKING








































EN	ER	PICTURE	REF	PRODUCT NAME	SPECIFICATION
IMPRESSION POST					
			320EN0A0	Impression Post EN	Aluminium Cap
			320ER0A1	Impression Post ER	Aluminium Cap
			690NA029	Impression Post Caps	Aluminium Kit 3pcs
			325EN0A0	Impression Post Pro EN	Plastic Cap
			325ER0A0	Impression Post Pro ER	Plastic Cap
			690NA091.10	Caps Kit Pro	Ø5.1mm Kit 10pcs
			690NA068	Impression Post Screw	M1.8 HEX1.20 H7.9mm
			690NA067	Impression Post Screw	M2 HEX1.20 H8.1mm

































EN	ER	PICTURE	REF	PRODUCT NAME	SPECIFICATION
IMPRESSION POST PROPICK-UP					
●			323EN0A0	Impression Post Propick-Up EN	HUseful22.8mm Long screw
			323EN0A1	Impression Post Propick-Up EN	HUseful17.8mm Short screw
			323EN0R0	Impression Post Propick-Up EN	HUseful22.8mm Rotating Long screw
			323EN0R1	Impression Post Propick-Up EN	HUseful17.8mm Rotating Short screw
●			690NA222	Impression Post Pick-Up Screw	M1.8 HEX1.20 H26mm
			690NA190	Impression Post Pick-Up Screw	M1.8 HEX1.20 H21mm
	●		323ER0A0	Impression Post Propick-Up ER	HUseful21.5mm Long screw
			323ER0A1	Impression Post Propick-Up ER	HUseful16.5mm Short screw
			323ER0R0	Impression Post Propick-Up ER	HUseful21.5mm Rotating Long screw
			323ER0R1	Impression Post Propick-Up ER	HUseful16.5mm Rotating Short screw
	●		690NA223	Impression Post Pick-Up Screw	M2 HEX1.20 H24.4mm
			690NA191	Impression Post Pick-Up Screw	M2 HEX1.20 H19.4mm
IMPLANT REPLICA					
●			301EN0A0	Implant Replica EN	
	●		301ER0A1	Implant Replica ER	

INTERIM RESTORATIONS































































EN	ER	PICTURE	REF	PRODUCT NAME	SPECIFICATION
TEMPORARY ABUTMENTS					
●			210EN2A0	Temporary Abutment EN	
			210EN2R0	Temporary Abutment EN	Rotating
●			215EN2A0	Temporary Abutment EN	Peek
●			690NA012	Retentive Screw EN	M1.8 HEX1.20
	●		210ER2A0	Temporary Abutment ER	
			210ER2A1	Temporary Abutment ER	Rotating
	●		215ER2A0	Temporary Abutment ER	Peek
	●		690NA013	Retentive Screw ER	M2 HEX1.20

CEMENT-RETAINED PROSTHESIS

EN	ER	PICTURE	REF	PRODUCT NAME	SPECIFICATION
STRAIGHT ABUTMENTS Ø 4.5 mm					
			220EN1A0	Straight Abutment EN	H1mm Ø4.5mm
			220EN2A0	Straight Abutment EN	H2mm Ø4.5mm
			220EN3A0	Straight Abutment EN	H3mm Ø4.5mm
			220EN4A0	Straight Abutment EN	H4mm Ø4.5mm
STRAIGHT ABUTMENTS Ø 5 mm					
			220ER2A1	Straight Abutment ER	H2mm Ø5mm
			220ER4A0	Straight Abutment ER	H4mm Ø5mm
STRAIGHT ABUTMENTS Ø 6 mm					
			220ER2A3	Straight Abutment ER	H2mm Ø6mm
			220ER4A1	Straight Abutment ER	H4mm Ø6mm
STRAIGHT ABUTMENTS Ø 7.5 mm					
			220ER2A2	Straight Abutment ER	H2mm Ø7.5mm
			220ER4A2	Straight Abutment ER	H4mm Ø7.5mm
ESTHETIC ABUTMENT					
			219EN0A0	Esthetic Abutment EN	
			219ER0A0	Esthetic Abutment ER	
ANGLED ABUTMENTS 15°					
			220EN2D0	Angled Abutment EN	15° H2mm
			220EN4D0	Angled Abutment EN	15° H4mm
			220ER2D0	Angled Abutment ER	15° H2mm
			220ER4D0	Angled Abutment ER	15° H4mm
ANGLED ABUTMENTS 25°					
			220EN2F0	Angled Abutment EN	25° H2mm
			220EN4F0	Angled Abutment EN	25° H4mm
			220ER2F0	Angled Abutment ER	25° H2mm
			220ER4F0	Angled Abutment ER	25° H4mm
RETENTIVE SCREWS FOR ABUTMENTS					
			690NA012	Retentive Screw	M1.8 HEX1.20
			690NA013	Retentive Screw	M2 HEX1.20
			VLE14TIT	Lingual Screw Esthetic Abut.	M1.4 HEX0.9

EN	ER	PICTURE	REF	PRODUCT NAME	SPECIFICATION
BT LINK					
●			246EN1A0	BT Link EN	H1mm Ø4.1mm
			246EN1A1	BT Link EN	H1mm Ø4.1mm Rotating
●			247EN1A0	Base BT Link EN	H1mm Ø4.1mm no Cap.
			247EN1A1	Base BT Link EN	H1mm Ø4.1mm no Cap. Rot.
●			205NA001.05	Castable Plastic Abut. BT Link	H1mm Ø4.7mm Kit 5pcs
●			690NA077	Retentive Screw EN	M1.8 HEX1.20 H7.8mm TP
	●		246ER1A0	BT Link ER	H1mm Ø4.6mm
			246ER1A1	BT Link ER	H1mm Ø4.6mm Rotating
	●		247ER1A0	Base BT Link ER	H1mm Ø4.6mm no Cap.
			247ER1A1	Base BT Link ER	H1mm Ø4.6mm no Cap. Rot.
	●		205NA002.05	Castable Plastic Abut. BT Link	H1mm Ø5.2mm Kit 5pcs
	●		244ER0A0	SIR Link ER	H0.8mm Ø4.1mm
	●		690NA076	Retentive Screw ER	M2 HEX1.20 H7.5mm TP
BT GRIP					
	●		530JD036	Screwdriver JD BT GRIP	HEX1.50 L15 mm (Short)
			530JD037	Screwdriver JD BT GRIP	HEX1.50 L30 mm (Long)
	●		248ER1A0	X3 Link ER	H1mm Ø4.6mm Multi High
	●		690NA236	Retentive Screw BT GRIP ER	M2 HEX 1.50
CAST-ON TECHNIQUE					
●			245EN1A0	Gold Abutment EN	H1mm
			245EN1R0	Gold Abutment EN	H1mm Rotating
●			240EN1A0	CoCr Abutment EN	H1.5mm
			240EN1R0	CoCr Abutment EN	H1.5mm Rotating
●			690NA012	Retentive Screw EN	M1.8 HEX1.20
●			690NA018	Retentive Screw Gold EN	M1.8 HEX1.20
	●		245ER1A0	Gold Abutment ER	H1mm
			245ER1R0	Gold Abutment ER	H1mm Rotating
	●		240ER1A0	CoCr Abutment ER	H1.5mm
			240ER1R0	CoCr Abutment ER	H1.5mm Rotating
	●		690NA013	Retentive Screw ER	M2 HEX1.20
	●		690NA019	Retentive Screw Gold ER	M2 HEX1.20
CASTABLE PLASTIC ABUTMENTS					
●			205EN1A1	Castable Plastic Abutment EN	
			205EN1A1.10	Castable Plastic Abutment EN	Kit 10pcs
●			205EN1R1	Castable Plastic Abutment EN	Rotating
			205EN1R1.10	Castable Plastic Abutment EN	Rotating Kit 10pcs
●			690NA012	Retentive Screw EN	M1.8 HEX1.20
●			690NA018	Retentive Screw Gold EN	M1.8 HEX1.20
	●		205ER2A0	Castable Plastic Abutment ER	
			205ER2A0.10	Castable Plastic Abutment ER	Kit 10pcs
	●		205ER2A1	Castable Plastic Abutment ER	Rotating
			205ER2A1.10	Castable Plastic Abutment ER	Rotating Kit 10pcs
	●		690NA013	Retentive Screw ER	M2 HEX1.20
	●		690NA019	Retentive Screw Gold ER	M2 HEX1.20
SCAN ABUTMENT					
●			351EN1A0	Scan Abutment Extra-oral EN	
			352EN1A0	Scan Abutment Intra-oral EN	
●			690NA077	Retentive Screw EN	M1.8 HEX1.20 H7.8mm TP
	●		351ER1A0	Scan Abutment Extra-oral ER	
			352ER1A0	Scan Abutment Intra-oral ER	
	●		690NA076	Retentive Screw ER	M2 HEX1.20 H7.5mm TP

SCREW-RETAINED PROSTHESIS

EN	ER	PICTURE	REF	PRODUCT NAME	SPECIFICATION
BT4 STRAIGHT ABUTMENTS					
			265EN1R0	BT4 Angled Abutments EN	Rotating H1mm
			265EN2R0	BT4 Angled Abutments EN	Rotating H2mm
			265EN3R0	BT4 Angled Abutments EN	Rotating H3mm
			265ER1R1	BT4 Angled Abutments ER	Rotating H1mm
			265ER2R1	BT4 Angled Abutments ER	Rotating H2mm
			265ER3R1	BT4 Angled Abutments ER	Rotating H3mm
BT4 ANGLED ABUTMENTS 17°					
			266EN2L0	BT4 Angled Abutments EN	17° H2mm Ø4.8mm
			266EN3L0	BT4 Angled Abutments EN	17° H3mm Ø4.8mm
			690NA066	Retentive Screw BT4	M1.8 Angled Abutments IR-EN
			266ER2L0	BT4 Angled Abutments ER	17° H2mm Ø4.8mm
			266ER3L0	BT4 Angled Abutments ER	17° H3mm Ø4.8mm
			690NA065	Retentive Screw BT4	M2 Angled Abutments ER
BT4 ANGLED ABUTMENTS 30°					
			266EN3G0	BT4 Angled Abutments EN	30° H3mm Ø4.8mm
			690NA066	Retentive Screw BT4	M1.8 Angled Abutments IR-EN
			266ER3G0	BT4 Angled Abutments ER	30° H3mm Ø4.8mm
			690NA065	Retentive Screw BT4	M2 Angled Abutments ER
BT4 CAPS					
			330NA0A0.04	Covering Caps BT4	H5 Kit 4pcs
			690NA024	Retentive Screw	M1.4 HEX1.20 10N
INTERIM RESTORATION BT4					
			267NA0A0	BT4 Titanium Abutment	
			207NA0A0	Castable Plastic Abutment BT4	
			207NA0A1	Castable Plastic Abutment BT4	No screw
			690NA024	Retentive Screw	M1.4 HEX1.20 10N
			311NA0A0	Impression Post Pick-up BT4	with Long Screw
			690NA031	Impression Post Pick-Up Screw	M1.4 HEX1.20 H17mm
			303NA0A0	Abutment Replica BT4	
BT4 SCAN ABUTMENT					
			351BT1A1	Scan Abutment Extra-oral BT	Rotating
			352BT1A1	Scan Abutment Intra-oral BT	Rotating
			690NA024	Retentive Screw	M1.4 HEX1.20 10N
BT4 BT LINK					
			246BT1A1	BT Link BT	H1mm Ø4.8mm Rotating
			247BT1A1	Base BT Link BT	H1mm Ø4.8mm no Cap. Rot.
			205NA003.05	Castable Plastic Abut. BT Link	H1mm Ø5.4mm Kit 5pcs
			690NA024	Retentive Screw	M1.4 HEX1.20 10N
CAST-ON TECHNIQUE					
			240BT1R0	CoCr Abutment BT	H1.5mm Rotating
			690NA024	Retentive Screw	M1.4 HEX1.20 10N

OVERDENTURE


SPHERO®



	SPHERO® Block Normo Sphere Ø 2.5mm		SPHERO® Block Micro Sphere Ø 1.8mm		SPHERO® Flex Sphere Ø 2.5mm 0°-7.5° Divergence	
	REF	SPECIFICATION	REF	SPECIFICATION	REF	SPECIFICATION
KR	254KR1A0	H1mm	255KR1A0	H1mm	256KR1A0	H1mm
	254KR2A0	H2mm	255KR2A0	H2mm	256KR2A0	H2mm
	254KR3A0	H3mm	255KR3A0	H3mm	256KR3A0	H3mm
	254KR4A0	H4mm	255KR4A0	H4mm	256KR4A0	H4mm
	254KR5A0	H5mm	255KR5A0	H5mm	256KR5A0	H5mm
	254KR6A0	H6mm	255KR6A0	H6mm	256KR6A0	H6mm
	254KR7A0	H7mm	255KR7A0	H7mm	256KR7A0	H7mm
IR	254IR1A0	H1mm	255IR1A0	H1mm	256IR1A0	H1mm
	254IR2A0	H2mm	255IR2A0	H2mm	256IR2A0	H2mm
	254IR3A0	H3mm	255IR3A0	H3mm	256IR3A0	H3mm
	254IR4A0	H4mm	255IR4A0	H4mm	256IR4A0	H4mm
	254IR5A0	H5mm	255IR5A0	H5mm	256IR5A0	H5mm
	254IR6A0	H6mm	255IR6A0	H6mm	256IR6A0	H6mm
	254IR7A0	H7mm	255IR7A0	H7mm	256IR7A0	H7mm
IW	254IW1A0	H1mm	255IW1A0	H1mm	256IW1A0	H1mm
	254IW2A0	H2mm	255IW2A0	H2mm	256IW2A0	H2mm
	254IW3A0	H3mm	255IW3A0	H3mm	256IW3A0	H3mm
	254IW4A0	H4mm	255IW4A0	H4mm	256IW4A0	H4mm
	254IW5A0	H5mm	255IW5A0	H5mm	256IW5A0	H5mm
	254IW6A0	H6mm	255IW6A0	H6mm	256IW6A0	H6mm
	254IW7A0	H7mm	255IW7A0	H7mm	256IW7A0	H7mm
EN	254EN1A0	H1.4mm	255EN1A0	H1mm	256EN1A0	H1mm
	254EN2A0	H2mm	255EN2A0	H2mm	256EN2A0	H2mm
	254EN3A0	H3mm	255EN3A0	H3mm	256EN3A0	H3mm
	254EN4A0	H4mm	255EN4A0	H4mm	256EN4A0	H4mm
	254EN5A0	H5mm	255EN5A0	H5mm	256EN5A0	H5mm
	254EN6A0	H6mm	255EN6A0	H6mm	256EN6A0	H6mm
	254EN7A0	H7mm	255EN7A0	H7mm	256EN7A0	H7mm
ER	254ER1A0	H1mm	255ER1A0	H1mm	256ER1A0	H1mm
	254ER2A0	H2mm	255ER2A0	H2mm	256ER2A0	H2mm
	254ER3A0	H3mm	255ER3A0	H3mm	256ER3A0	H3mm
	254ER4A0	H4mm	255ER4A0	H4mm	256ER4A0	H4mm
	254ER5A0	H5mm	255ER5A0	H5mm	256ER5A0	H5mm
	254ER6A0	H6mm	255ER6A0	H6mm	256ER6A0	H6mm
	254ER7A0	H7mm	255ER7A0	H7mm	256ER7A0	H7mm

Accessori Sphero®

PICTURE	REF	PRODUCT NAME	SPECIFICATION
	530JD030	Wrench Cric. Sphero Flex/Block	Rhein83® 771CEF

NOTA Every SPHERO® as listed above includes the following products: 1pc. Titanium Abutment with self-aligning 2.5mm or 1.8 mm sphere, 2pcs. Soft Retention Pink Caps, 1pc. Stainless Steel Housing, 1pc. Protective Disk and 3 pcs. Directional Rings.
These Devices are produced by Rhein83.  s.r.l. Via E. Zago, 10/ABC, 40128 Bologna Italy.



Abutment Connector

KR		IR		ER		EN	
REF	SPECIFICATION	REF	SPECIFICATION	REF	SPECIFICATION	REF	SPECIFICATION
261KR1A0	H1mm	261IR1A0	H1mm	261ER1A0	H1mm	261EN1A0	H1mm
261KR2A0	H2mm	261IR2A0	H2mm	261ER2A0	H2mm	261EN2A0	H2mm
261KR3A0	H3mm	261IR3A0	H3mm	261ER3A0	H3mm	261EN3A0	H3mm
261KR4A0	H4mm	261IR4A0	H4mm	261ER4A0	H4mm	261EN4A0	H4mm
261KR5A0	H5mm	261IR5A0	H5mm	261ER5A0	H5mm	261EN5A0	H5mm
261KR6A0	H6mm	261IR6A0	H6mm	261ER6A0	H6mm	261EN6A0	H6mm

Connector and Locator® Accessories

PICTURE	REF	PRODUCT NAME	SPECIFICATION
	530JD029	Screwdriver JD Connector	L10mm (for reversible torque wrench JD)
	530HS027	Handpiece Driver Connector	L 23mm
	530HS028	Handpiece Driver Connector	L 29mm
	540MA026	Angle Measur. Guide Connector	
	540MA027.04	Parallel Pin for Connector	Kit 4pcs
	321NA0A0	Transfer Locator®	
	301C00A0	Implant Replica Connector Abut.	
	690NA011	Replacement Males Locator®	
	690NA022	Kit Locator®	Metal Cap + Blockout Spacer
	690NA054.04	Replacement Males Locator®	0N Black Kit 4pcs
	690NA006.04	Replacement Males Locator®	15N Blue Kit 4pcs (700gr. 0-20°)
	690NA008.04	Replacement Males Locator®	30N Pink Kit 4pcs (1400gr. 0-20°)
	690NA010.04	Replacement Males Locator®	50N Neutral Kit 4pcs (2300gr. 0-20°)
	690NA005.04	Replacement Males Locator®	10N Red Kit 4pcs (700gr. 20-40°)
	690NA007.04	Replacement Males Locator®	20N Orange Kit 4pcs (900gr. 20-40°)
	690NA009.04	Replacement Males Locator®	40N Green Kit 4pcs (1400gr. 20-40°)
	690NA134.04	Replacement Males Locator®	0N Gray Kit 4pcs (0gr.)
	502MA024	Connector Core Tool 3 in 1	
	502MA025	Connector Male Removal Tip End	
	690NA320	Retaining Sleeve Connector	

NOTA Every "Connector Abutment" as listed above includes the following products: 1 pc. Connector Abutment; 1 pc. Denture Male Cap (Housing); 1 pc. Block-Out Spacer, 1 pc. each LOCATOR® Replacement Males (blue / pink / clear). LOCATOR® is a registered trademark of Zest Anchors Inc.

MATERIAL SPECIFICATIONS

TITANIUM GRADE 4 IMPLANTS

CHEMICAL COMPOSITION:	MAXIMUM VALUES (%)	TOLERANCE
Nitrogen (N)	0.05	+/- 0.02
Carbon (C)	0.08	+/- 0.02
Hydrogen (H)	0.015	+/- 0.002
Iron (Fe)	0.50	+/- 0.10 (%<0.25) +/- 0.15 (%>0.25)
Oxygen (O)	0.40	+/- 0.02 (%<0.20) +/- 0.03 (%>0.20)
Titanium (Ti)	balance	-

MECHANICAL PROPERTIES:	MINIMUM VALUES
Tensile stress:	550 MPa
Yield strength (0.2%):	483 MPa
Elongation at yield:	15 %
Section reduction:	25 %

This technical information complies with the express specification of the regulations in force for the use of grade 4 titanium in implantology:

- ASTM F67: Standard Specification for unalloyed titanium, for surgical implant applications.
- ISO 5832-2: Implant for surgery – Metallic Materials – Part 2: Unalloyed titanium.

TITANIUM GRADE 5 PROSTHETICS AND MINI IMPLANTS

CHEMICAL COMPOSITION:	MAXIMUM VALUES (%)	TOLERANCE
Nitrogen (N)	0.05	+/- 0.02
Carbon (C)	0.08	+/- 0.02
Hydrogen (H)	0.012	+/- 0.002
Iron (Fe)	0.25	+/- 0.10
Oxygen (O)	0.13	+/- 0.02
Aluminium (Al)	5.50-6.50	+/- 0.40
Vanadium (V)	3.50-4.50	+/- 0.15
Titanium (Ti)	balance	-

MECHANICAL PROPERTIES:	MINIMUM VALUES
Tensile stress:	860 MPa
Yield strength (0.2%):	795 MPa
Elongation at yield:	10 %
Section reduction:	25 %

This technical information complies with the express specification of the regulations in force for the use of grade 5 titanium in implantology:

- ASTM F136: Standard Specification for wrought Titanium-6Aluminium-4Vanadium ELI (Extra low Interstitial) Alloy for surgical implant applications;
- ISO 5832-3: Implant for surgery – Metallic Materials – Part 3: Wrought titanium 6-aluminium 4-vanadium alloy.

COBALT CHROME COBALT CCM®

CHEMICAL COMPOSITION: (%)	
Carbon (C)	max. 0.14
Silicon (Si)	max. 1.00
Manganese (Mn)	max. 1.00
Chromium (Cr)	26.00-30.00
Molybdenum (Mo)	5.00-7.00
Nickel (Ni)	max. 1.0
Iron (Fe)	max. 0.75
Nitrogen (N)	max. 0.25
Cobalt (Co)	balance

MATERIAL NO. AND NORMS	
DIN	CoCr28Mo
ISO	5832-12
AFNOR	CoCr28Mo
ASTM	F1537 alloy 1
UNS	R31537

MECHANICAL PROPERTIES	
Coefficient of Expansion (CTE)	13.2•10 ⁻⁶ °C ⁻¹
Melting range	1340-1440°C
Yield strength (R0.2)	up to 1115 MPa
Young Modulus E	241 GPa
Hardness	up to 46 HRC

PRECIOUS ALLOY FOR ABUTMENTS

COMPOSITION:	
Gold (Au)	60.0 %
Platinum (Pt)	24.9 %
Palladium (Pd)	15.0 %
Iridium (Ir)	0.1 %
PHYSICAL AND MECHANICAL PROPERTIES:	
Density:	18.1 g/cm ³
Melting range:	1350 – 1460 °C
Coefficient of Expansion (CTE) 25-500°C – 25-600°C:	12.7•10 ⁻⁶ °C ⁻¹ – 12.9•10 ⁻⁶ °C ⁻¹
Modulus of elasticity (tensile test):	110 GPa
Elongation at yield:	18 – 12 %
Breaking load:	580 – 810 MPa
Yield strength (0.2%):	450 – 720 MPa
Vickers Hardness HV5/30:	150 – 205 – 230

PRECIOUS ALLOY FOR GOLD RETENTIVE SCREWS

COMPOSITION:	MAXIMUM VALUES (%)	TOLERANCE
Gold (Au)	0,5	+/- 0.2
Gallium (Ga)	2	+/- 0.2
Copper (Cu)	10	+/- 0.5
Iridium (Ir)	7	+/- 0.5
Ruthenium (Ru)	0.03	+/- 0.02
Rutenio (Ru)	0.1	+/- 0.09
Palladium (Pd)	balance	

MECHANICAL PROPERTIES:	MINIMUM VALUES (%)
Tensile stress:	586 - 862 MPa
Yield strength (0.2%):	483 - 690 MPa
Elongation:	5 - 20 %
Young's Modulus:	138 GPa

PHYSICAL PROPERTIES:	
Melting Range	1450 – 1500 °C
Coefficient of Expansion (CTE) 25-500°C – 25-600°C:	12.3•10 ⁻⁶ °C ⁻¹

The temporary abutments in PEEK and the SCAN ABUTMENT are made of PEEK / TECAPEEK CLASSIC (chemical name Polietereeterketone). This material is suitable to stay in contact with tissue for up to 180 days.

Depending on the intended use, the Biotec instrumental is made of specific types of stainless steel.

SYMBOLS USED ON LABELS



Legal manufacturer



Use-by date: indicates the date after which this device is not to be used



Products with the CE mark in accordance with Directive 93/42/EEC and following modifications/integrations



Do not use if packaging is damaged



0426
Number of the notification body



Do not reuse



Consult instructions for use



Keep away from sunlight



ifu.btk.dental

Electronic instructions for use available online ifu.btk.dental



Sterile by gamma irradiation



Caution; see instructions for use



Catalogue number



Lot/batch number

DELIVERY TERMS & CONDITIONS

RESPONSABILITY

The use of BTK medical devices is reserved exclusively for personnel with the necessary qualifications for the exercise. An improper or incorrect use of the devices can cause the failure or worse, injury to the patient or the user. BTK implant systems should only be used with original BTK components and instruments and in accordance with the specific BTK instructions. Combining with different devices may cause a failure. Biotec must not and can not control the procedures for using the product for implant-prosthetic treatment. Therefore, Biotec assumes no responsibility for the application of the device and its processing nor for any incongruous use of the device under the surgical or prosthetic profile, nor in any case for failure, adverse reactions or damage to the patient or dentist as a result of application of the product.

STERILITY OF WARRANTY AND DISPOSABLE

Dental implants are supplied STERILE (gamma ray sterilization). The sterility of the medical implant is guaranteed only according to the following conditions: the expiry date stated on the packaging is still valid; there is a red dot on the sterile vial which demonstrates that it has undergone gamma ray irradiation; the sealed package has not been opened and does not show any signs of damage. Compliance with all these conditions must be ensured; alternatively do not use the device.

Surgical components, laboratory accessories and instruments are not supplied in sterile packs, therefore before use they must be properly CLEANED and STERILIZED, as shown in the instructions for use. Biotec dental implants, prosthetics and laboratory accessories are designed for SINGLE USE. In fact, reuse is a potential risk and could damage the construction of the device, making it inappropriate for its intended use. Biotec explicitly declares the single-use of MD and assumes no responsibility for any re-use by users.

STORAGE

Biotec products must be stored at room temperature and protected from direct heat or sunlight and dust.

INSTRUCTIONS FOR USE

The information in this manual is not intended to be exhaustive for BTK implant systems. It is recommended that new customers follow the training courses that Biotec makes available with trained personnel and clinicians who are experts in implantology and in the use of BTK devices. The complete and updated user manuals, which allow the correct use of the product, are available online (www.btk.dental) or at BTK and / or the local distributor.

AVAILABILITY

Not all products described here are available in ExtraEU countries. For more information, please contact BTK and / or your local distributor.

RETURNS

Biotec does not accept returned goods if the packaging seals are broken or not conforming to the sale specifications of the company.

GUARANTEE

We constantly guarantee that the quality of our products and services meets the high expectations of our customers and their patients. Specialized professionals are committed to offering complete solutions in applied research, engineering, training and related activities. Biotec is available to customers in the event that a defect in the product or its use is found.

VALIDITY

The contents are updated at the date of publication. This manual replaces all previous editions.

CASE DOCUMENTATION AND TRACEABILITY

BTK absolutely recommends documenting implant cases comprehensively at the clinical, radiographic, photographic and statistical levels. The clinician must guarantee the traceability of the devices used. It is advisable to use the adhesive labels included in the packaging of the BTK devices, which show the code and lot of the device used, for the purpose of documentation on the medical records and on the relative implant passport of the patient.

TRAINING

Comprehensive and regular training ensures long-term implant success.

Be advised that we strongly recommend regular education events in order to update your know-how and clinical expertise.

DELIVERY TERMS

BTK delivery terms are 1 working day for order received before 12.00 p.m. of the previous day in Italy; except for islands where delivery is evaluated to be 2 working days. For export deliveries contact Biotec offices.

QUALITY STANDARD

Owing to extensive research, development and to a strict quality standard, we guarantee premium quality materials and products. Our products meet the requirements of directive 93/42 /EEC and subsequent amendments and additions, and therefore have the CE mark, in accordance with the corresponding law. BTK has a quality system certified UNI EN ISO 9001 and UNI EN ISO 13485.

CAUTION

In addition to the instructions for use, warnings and risks reported both in this document and in the instructions for use, it must always be ensured that the devices used in the oral cavity are not aspirated or swallowed by the patient.

COPYRIGHT AND TRADEMARKS

BTK documents cannot be reprinted or published, in whole or in part, without the written permission of BTK.

BTK, BIOTEC, The Smile System® or Implanting Trust, Smile Again® are trademarks or registered trademarks of BIOTEC s.r.l. and / or its affiliates.



BTK PERSONAL TUTOR

A program for individual case planning and execution supported by experienced professionals in order to leverage know-how and maximize clinical experience with the aim to achieve sustainable high patient satisfaction rates.

BTK is always at your disposal for any request for further follow-up or information, promoting periodic and ad-hoc training course.

CERTIFIED QUALITY SYSTEM

BIOTEC is certified UNI EN ISO 9001
and UNI EN ISO 13485.



CE marked product, in accordance with Directive 93/42/EEC and subsequent modifications and additions.

MADE IN ITALY USED GLOBALLY



We constantly ensure that the quality of our products and services meet the high expectations of our customers and their patients.

Specialized professionals are taking care to offer comprehensive solutions in applied research, engineering, education and related activities.