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**i** The following information supplements the CS 1 Installation and Operating Instructions.

Always refer to the Installation and Operating Instructions of the unit in addition, as these contain important information concerning safety instructions, set-up, electrical connections as well as disinfection and cleaning etc..

**1. Correct Usage**

The CS 1 with heavy particle collector has been designed to be integrated into treatment stations where no amalgam will be used or encountered. The main task of the heavy particle collector is to collect the larger, heavier soluble particles, e.g. calcium carbonate as used for teeth cleaning. This is designed to prevent this powder from getting into the waste system and the plumbing line where it would first be deposited as a sediment and in the long term would lead to blockage within the system.

**2. Function**

The heavy particle collector works according to the sedimentation principle. Fluids drain from the spittoon of the CS1 (6) and then flow through the inlet hose (1) into the collection vessel (2). The max. fluid flow rate must not exceed 3 l/min. A filter cartridge within the collector vessel first brings the fluid to rest. The special design of the housing prevents the fluids from proceeding directly into the waste water connection (3), and ensures instead that there is a longer slowing and sedimentation period.

- 4 Air extraction
- 5 Suction connection
- 7 Emergency drain

**3. Maintenance**

**In order to prevent any danger of infection we strongly recommend that during the collector vessel changing water-proof protective gloves are worn.**

Maintenance of the heavy particle collector is reduced to regular inspection, at least once a week, with either changing or cleaning of the collector. When the collection vessel is cleaned then the filter insert should be replaced by a new one. General wear and tear of the O-ring collector seal can lead to leaks, so this should also be checked and replaced when necessary.

**Disposable materials:**

- Collection vessel with O-ring, 5 filter cartridges 7117-034-00
- O-Ring for collection vessel, 10 filter cartridges 7117-035-00

**4. Fault**

**In order to prevent any danger of infection we strongly recommend that during the collector vessel changing water-proof protective gloves are worn.**

If the collection vessel is not regularly checked or changed, this can lead to the collector vessel becoming too full causing a blockage in the system. In this case loud gurgling noises will be heard at the CS 1 or there will be a back flow of fluids into the spittoon. Should this arise the emergency drainage hose can be used (7) to drain the system to such an extent that the collector vessel can then be removed safely with no danger of spillage or overflow.