

General Instructions

Range of table curing lamps

This document is an English translation of the original French version.
Reference J02900 version V7 and drawing number RG41FR040G

Contents

| | |
|---|-----------|
| 1 Documentation | 1 |
| 1.1 Associated documentation | 1 |
| 1.2 Electronic documentation | 1 |
| 2 Warnings | 3 |
| 2.1 Federal Law | 3 |
| 2.2 Warning applicable to all countries in which the device is sold | 3 |
| 2.3 User population | 3 |
| 2.4 Specific user training | 3 |
| 2.5 Patient population | 3 |
| 2.6 Patient population restriction | 3 |
| 2.7 Parts of the body or types of tissues treated | 3 |
| 2.8 Applied parts | 3 |
| 2.9 Essential performance | 3 |
| 2.10 Basic safety in normal use | 4 |
| 2.11 Normal usage conditions | 4 |
| 2.12 Service life | 4 |
| 3 Interactions, contraindications, prohibitions | 5 |
| 3.1 Interference with other medical devices | 5 |
| 3.2 Using accessories not supplied by the manufacturer | 5 |
| 3.3 Prohibited uses | 5 |
| 3.4 Moving the medical device | 5 |
| 3.5 Precautions for use | 5 |
| 4 Electromagnetic compatibility | 7 |
| 4.1 Cable length | 7 |
| 4.2 Recommended separation distances | 7 |
| 4.3 Electromagnetic emissions | 7 |
| 4.4 Magnetic and electromagnetic immunity | 8 |
| 4.5 Electromagnetic immunity, handheld radiofrequency equipment | 9 |
| 5 Disinfection and sterilising | 11 |
| 6 Regulations and standards | 13 |
| 6.1 Latest document update | 13 |
| 6.2 Manufacturer identification | 13 |
| 6.3 Manufacturer responsibility | 13 |
| 6.4 Branch addresses | 13 |
| 6.5 Disposal and recycling | 15 |
| 6.6 Applicable standards and regulations | 15 |
| 6.7 Symbols | 15 |
| 7 Glossary | 17 |
| 8 Index | 19 |

1 Documentation

This document contains the following information:

- Patient, practitioner and environment safety
- Installing your medical device in optimum conditions
- Identifying the manufacturer or the latter's representatives if necessary

1.1 Associated documentation

| Document title | References |
|---|------------|
| Cleaning, disinfection and sterilisation instructions for the MINILED optical guide | J02941 |
| Cleaning and disinfection instructions for the MINILED rigid protection shield | J05541 |
| Cleaning and disinfection instructions for the MINILED flexible protection shield | J05551 |
| Consulting electronic user instructions | J00007 |
| MINILED Blue Ray 3 User manual | J05251 |
| MINILED ACTIVE User manual | J05271 |
| MINILED STANDARD User Manual | J02541 |
| MINILED SUPERCHARGED User manual | J02261EN |
| MINILED ORTHO 2 User manual | J05221 |

1.2 Electronic documentation



The user instructions for your device are available in electronic format on the specified website and not in printed format. If the website is unavailable, try again later. You can also request a free printed copy of the user instructions within seven days via our website, by telephone or in writing.

The electronic user instructions are available in PDF format (Portable Document Format). You will need to have a PDF file reader installed to read the electronic user instructions. It is important for you to have read and understood the content of the user instructions relating to the use of your device and its accessories.

| Never use your device without first reading the user instructions.

The device user instructions can be consulted at www.satelec.com/documents

When you receive your device, you are asked to print and/or to download all documents or sections of documents that you may need to refer to in the event of an emergency, if you are unable to connect to the internet or if your electronic display tool is not working (computer, tablet, etc.). We recommend that you visit the website regularly to consult and to download the latest version of your device's user instructions. Users are asked to keep documentation close at hand for reference when necessary.

All paper or electronic documentation relating to your medical device must be kept for the device's entire service life. Keep the original documentation for your medical device and its accessories for future reference. When loaning out or selling the medical device, the documentation must be provided with it.

2 Warnings

2.1 Federal Law

| The indication below applies to the United States of America only.

The United States Federal Law restricts the use of this medical device in its territory to qualified, fit and certified dental health professionals (either directly or under their supervision).

2.2 Warning applicable to all countries in which the device is sold

| The information below is based on the requirements of standards to which the manufacturers of medical devices must adhere (as stated in standard IEC62366).

2.3 User population

Users must know and comply with the rules of dental practice in compliance with knowledge acquired in the field and the key medical hygiene principles including cleaning, disinfection and sterilisation of medical devices.

The medical device can be used by any adult dental practitioner of any weight, age, height, gender or nationality.

The user must wear gloves.

The user is not the patient.

The user must not be prone to any of the following:

- visual impairments: any vision problems must be corrected by glasses or lenses.
- arm disability that may prevent the user from holding a handpiece;
- hearing difficulties that could prevent the user hearing audible alarms depending on medical devices;
- difficulty memorizing or concentrating that could affect the setting of sequences or the performance of treatment protocols.

2.4 Specific user training

No specific training other than initial professional training is required to use this medical device.

The practitioner is responsible for performing clinical treatments and for dangers that may arise due to a lack of skill and/or training.

2.5 Patient population

This medical device is designed to be used with the following patient populations:

- children;
- Teenagers,
- Adults,
- Old Age Pensioners.

This medical device can be used on any patient of any weight (except children), age, height, gender or nationality.

2.6 Patient population restriction

This medical device must not be used on the following patient populations:

- infants ;
- pregnant or breastfeeding women due to restrictions associated with the possible use of medical solutions such as anaesthetics;
- patients with medical issues,
- Patients allergic to some of the medical device components,
- patients with a clinical site not suitable for treatment,

The patient must be calm, relaxed, still, ideally lying flat on a dental chair.

The user is the only person who can decide whether or not to treat his/her patients.

2.7 Parts of the body or types of tissues treated

Treatment must only be performed on the patient's oral environment.

2.8 Applied parts

| | |
|---|---------------|
| Part in indirect contact with the patient | Optical guide |
|---|---------------|

2.9 Essential performance

As stated in the applicable safety standard pertaining to electrical medical devices, The manufacturer has determined that the medical device did not manage essential performances.

2.10 Basic safety in normal use

The active part, the handpiece is held by the practitioner throughout the treatment.

As a highly skilled medical expert, the practitioner can immediately detect any problems at the treatment area and react accordingly.

It is advisable to have a spare medical device or an alternative means with which to perform the medical treatment in the event of device failure.

2.11 Normal usage conditions

The normal usage conditions are as follows:

- storage;
- installation;
- use;
- maintenance;
- disposal.

2.12 Service life

Systematically replace an optical guide that is damaged due to accidental impact. Do not use the MINILED if it or the optical guide are damaged or faulty.

The control elements are sensitive to pressure and to wear that may be caused by excessive pressure.

3 Interactions, contraindications, prohibitions

This includes information relating to the interactions, contraindications and prohibited operations known by the manufacturer on the date on which this document was written.

3.1 Interference with other medical devices

The medical device presents potential risks due to the emission of electromagnetic fields. Interferences may occur when the system is used on patients fitted with implantable medical devices such as a pacemaker, deep brain stimulator or vagus nerve stimulator.



It may in particular cause malfunction of all types of active implanted device:

- before using this medical device, check whether patients and practitioners are fitted with a device of this type (active or inactive);
- explain the situation;
- weigh up the benefits versus the risks and contact your patient's cardiologist or another qualified health professional prior to starting treatment;
- keep this system away from implantable devices;
- apply suitable emergency measures and act fast if the patient shows signs of being unwell.

Symptoms such as an increased heart beat, irregular pulse or dizziness may indicate a malfunction of a pacemaker or an implantable defibrillator.

The medical device is not designed to withstand electrical defibrillation shocks.

3.2 Using accessories not supplied by the manufacturer

The medical device was designed and developed with its accessories to guarantee maximum safety and performance. The use of accessories from another source could put you and your patients at risk and could damage your medical device.

Do not try to connect accessories not provided by SATELEC, a company of Acteon group to your medical device connector (s) or to the handpiece.

Even if the manufacturer or dealer of your accessory claims full compatibility with SATELEC, a company of Acteon group equipment, it is advisable to exercise caution with regards to the origin and safety of the product offered. Look out in particular for lack of information, information in a foreign language, very attractive prices, suspect appearance, mediocre quality or premature wear. If necessary, contact an approved dealer or the SATELEC, a company of Acteon group Customer Service team.

3.3 Prohibited uses

- Never cover the medical device and/or obstruct the air inlets.
- Do not immerse or use outdoors.
- Do not place the medical device next to a source of heat or in direct sunlight.
- Do not expose the medical device to water spray or mist.
- Do not use the medical device in an AP or APG gas-filled atmosphere.

The medical device is not designed to operate near a source of ionising radiation.

A hot/cold temperature contrast can cause condensation to form in the medical device, which may be dangerous. If the medical device needs to be moved from a cold place to a warm place, do not use the device immediately. Wait until it reaches room temperature.

The medical device may not be stored or used outside the temperature, atmospheric pressure and humidity ranges recommended in the User Manual supplied with your medical device.

Only use the medical device for the purpose for which it has been designed.

3.4 Moving the medical device

After its initial installation, the medical device is not designed to be moved. The medical device must be fixed to ensure that it cannot be removed or moved without the use of a tool.

3.5 Precautions for use

Do not use the MINILED on patients currently afflicted by or with a past history of the following conditions:

- photo-biological reactions, including solar urticaria or erythropoietic Protoporphyrria;
- treatment with photosensitizing medicines, including Methoxsalen and Chlorotetracycline.

Practitioners and patients who have suffered from retina or lens problems or who have had an operation on their eyes, in particular to treat cataracts, must consult their ophthalmologist before using the MINILED. Even with your ophthalmologist's approval, proceed with caution. The intensity of the light could cause accidents.

Throughout the duration of the procedure, the practitioner and the patient must wear class 2 safety goggles designed for use with medical devices emitting radiation with a wavelength of less than 500 nm.

Light rays emitted by the MINILED can be harmful and must never be aimed directly at the eyes even if the practitioner and the patient are wearing class 2 safety goggles.

Over exposure of the pulp and soft tissues to light rays may result in heat generation and may harm the patient. To prevent any feeling of heat, carry out the cure cycles with a rest time of 30 seconds between each 10-second cure cycle.

Do not use accessories that are damaged or that have unrecognisable markings, corrosion or sharp or cutting surfaces.

4 Electromagnetic compatibility

All the information below is based on the requirements of standards to which the manufacturers of electrical medical devices must adhere (as stated in standard IEC60601-1-2).

The medical device complies with the electromagnetic compatibility standards in force. However, the user must make sure that any electromagnetic interference does not create an additional risk, such as those created by radiofrequency transmitters, or other electronic devices.

This chapter contains the information required for you to install and use your medical device in optimum conditions in terms of electromagnetic compatibility.

The different medical device cords must be kept away from each other.

Some types of mobile telecommunication devices such as mobile phones may interfere with the medical device. The separation distances recommended in this chapter **MUST** be complied with.

The medical device must not be used near another device or placed on top of it. If this cannot be avoided, correct operation of the device in operating conditions must be checked prior to use.

The use of accessories other than those specified or sold by SATELEC, a company of Acteon group as replacement parts, may increase the transmission or reduce the immunity of the medical device.

4.1 Cable length

| Cables and accessories | Maximum length | Test type | In compliance with: |
|------------------------|----------------|--|---------------------|
| Cables/Cords | < 3 m | RF emission | CISPR 11, Class B |
| | | Harmonic current emission | IEC61000-3-2 |
| | | Voltage fluctuation and flickers | IEC61000-3-3 |
| | | Electrostatic discharge immunity | IEC61000-4-2 |
| | | Radiated immunity – Electromagnetic fields | IEC61000-4-3 |
| | | Electrical fast transient/burst immunity | IEC61000-4-4 |
| | | Surge immunity | IEC61000-4-5 |
| | | Immunity to conducted disturbances, induced by radiofrequency fields | IEC61000-4-6 |
| | | Radiated immunity - Magnetic fields | IEC61000-4-8 |
| | | Voltage dips, short interruptions and voltage variation immunity | IEC 61000-4-11 |

4.2 Recommended separation distances

The medical device is designed to be used in an electromagnetic environment in which interferences caused by radiofrequency radiation are controlled.

Do not use handheld radiofrequency communication devices within 30 cm (12 inches) of any part of the medical device, including its cables.

Aerial cables and external aeriels of handheld radiofrequency communication devices must not be positioned or used within 30 cm (12 inches) of any part of the medical device.

If the minimum distance is not adhered to when using handheld radiofrequency communication devices, this may impact the performance of the medical device.

4.3 Electromagnetic emissions

The medical device is designed for use in the electromagnetic environment described in the table below. The user and/or installer must ensure that the medical device is used in the environment described below.

| Emission test | Conformity | Electromagnetic environment - comments |
|--|------------|---|
| Electromagnetic radiation disturbance, radiated emissions (CISPR 11) | Group 1 | The medical device uses radiofrequency energy for its internal operation. Consequently, its radiofrequency emissions are very low and are not likely to create any interference with other nearby equipment. |
| Radiofrequency emission (CISPR 11) | Class A | The emission characteristics of this medical device make it suitable for use in industrial and hospital areas [Class A defined in CISPR 11]. When used in a residential environment, for which class B defined in CISPR 11 is normally required, this medical device may not provide adequate protection for radio frequency communications services. The user may need to take corrective measures such as the re-installation or reorientation of the medical device. |
| Radiofrequency emission (CISPR 11) | Class B | The medical device is suitable for use in a home-based health care setting or in a professional health care setting. |
| Harmonic current emission (IEC61000-3-2) | Class A | |
| Voltage fluctuation and flickers (IEC61000-3-3) | Compliant | |

The following medical devices are categorised as Class B radio-frequency equipment according to CISPR 11:

- MINILED ACTIVE
- MINILED STANDARD
- MINILED ORTHO 2
- MINILED SUPERCHARGED

4.4 Magnetic and electromagnetic immunity

The medical device is designed for use in the magnetic and electromagnetic environment described in the table below. The user and/or installer must ensure conformity of the electromagnetic environment.

| Immunity test | Test level in accordance with IEC60601 | Conformity level | Electromagnetic environment / comments |
|--|---|--|--|
| Electrostatic discharge (ESD) (IEC61000-4-2) | ± 8 kV on contact ± 15 kV in the air | ± 8 kV on contact ± 15 kV in the air | The medical device is suitable for use in a home-based health care setting or in a professional health care setting. |
| Electrical fast transient/burst immunity (IEC61000-4-4) | ± 2 kV for electricity supply lines ± 1 kV for signal ports Valid for medical devices with signal ports | ± 2 kV for electricity supply lines ± 1 kV for signal ports | The medical device is suitable for use in a home-based health care setting or in a professional health care setting. |
| Surge (IEC61000-4-5) | ±0.5 kV, ±1 kV between phases ±0.5 kV, ±1 kV, ±2 kV between phase and earth Valid for earthed medical devices | ±0.5 kV, ±1 kV between phases ±0.5 kV, ±1 kV, ±2 kV between phase and earth | The medical device is suitable for use in a home-based health care setting or in a professional health care setting. |
| Magnetic field at the assigned industrial frequency (IEC61000-4-8) | 30 A/m | 30 A/m | The magnetic field intensity must be equal to the level found in a home-based health care setting and in a professional health care establishment setting. |

| Immunity test | Test level in accordance with IEC60601 | Conformity level | Electromagnetic environment / comments |
|---------------------------------------|---|---|---|
| Voltage dip (IEC 61000-4-11) | 0% UT for 0.5 cycles at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% UT for 1 cycle and 70% UT for 25 cycles at 50 Hz for 30 cycles at 60 Hz Single phase at 0° | 0% UT for 0.5 cycles at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% UT for 1 cycle and 70% UT for 25 cycles at 50 Hz for 30 cycles at 60 Hz Single phase at 0° | The quality of the network supply must be equal to that of a home-based health care setting and a professional health care establishment setting. If the use of the system requires continuous operation during mains power outages, it is advisable to supply the medical device using a separate current source (UPS, etc.). |
| Voltage interruptions (IEC61000-4-11) | 0% UT for 250 cycles at 50 Hz for 300 cycles at 60 Hz | 0% UT for 250 cycles at 50 Hz for 300 cycles at 60 Hz | The quality of the network supply must be equal to that of a home-based health care setting and a professional health care establishment setting. If the use of the system requires continuous operation during mains power outages, it is advisable to supply the medical device using a separate current source (UPS, etc.). |

4.5 Electromagnetic immunity, handheld radiofrequency equipment

The medical device is designed for use in the magnetic and electromagnetic environment described in the table below. The user and/or installer must ensure conformity of the electromagnetic environment.

| Immunity test | Test level | Conformity level | Electromagnetic environment - comments |
|---|--|--|--|
| <p>Do not use handheld radiofrequency communication devices within 30 cm (12 inches) of any part of the medical device, including its cables.</p> <p>Aerial cables and external aeriels of handheld radiofrequency communication devices must not be positioned or used within 30 cm (12 inches) of any part of the medical device.</p> <p>If the minimum distance is not adhered to when using handheld radiofrequency communication devices, this may impact the performance of the medical device.</p> | | | |
| Radiated, radiofrequency, electromagnetic fields (IEC61000-4-3) | 10 V/m 80 MHz to 2.7 GHz 80% MA at 1 kHz | 10 V/m 80 MHz to 2.7 GHz 80% MA at 1 kHz | The medical device is suitable for use in a home-based health care setting or in a professional health care setting. |

| Immunity test | Test level | Conformity level | Electromagnetic environment - comments |
|---|---|---|--|
| Proximity fields transmitted by wireless radiofrequency communication devices (IEC 61000-4-3, temporary method) | 9 V/m 710 MHz, 745 MHz, 780 MHz, 5 240 MHz, 5 550 MHz, 5 785 MHz 27 V/m 385 MHz 28 V/m 450 MHz, 810 MHz, 870 MHz, 930MHz, 1 720 MHz, 1 845 MHz, 1 970 MHz, 2 450 MHz | 9 V/m 710 MHz, 745 MHz, 780 MHz, 5 240 MHz, 5 550 MHz, 5 785 MHz 27 V/m 385 MHz 28 V/m 450 MHz, 810 MHz, 870 MHz, 930MHz, 1 720 MHz, 1 845 MHz, 1 970 MHz, 2 450 MHz | The medical device is suitable for use in a home-based health care setting or in a professional health care setting. |
| Conducted disturbances, induced by radiofrequency fields (IEC61000-4-6) | 3 V/m 0.15 MHz to 80 MHz 6 V in the ISM band and bands between 0.15 MHz and 80 MHz, amateur radio bands included 80% MA at 1 kHz | 3 V/m 0.15 MHz to 80 MHz 6 V in the ISM band and bands between 0.15 MHz and 80 MHz, amateur radio bands included 80% MA at 1 kHz | The medical device is suitable for use in a home-based health care setting or in a professional health care setting. |

These specifications may not be applicable in all situations. The electromagnetic propagation is affected by the absorption and reflection of structures, objects and people.

The electromagnetic field intensity of fixed radiofrequency transmitters, such as base stations for portable phones (mobiles / wireless), mobile radios, radio amateurs, AM/FM radio transmissions and TV transmissions cannot be determined accurately by the theory.

To assess the electromagnetic environment caused by fixed radiofrequency transmitters, an electromagnetic environment measurement must be taken. If the measured intensity of the radiofrequency field in the product's immediate use environment exceeds the radiofrequency conformity level specified above, it is necessary to test product performance to check this complies with specifications. If abnormal performance is observed, additional measures may be necessary, such as changing the direction of or moving the product.

In the 150 kHz to 80 MHz frequency range, the electromagnetic fields must be less than 3 V/m.

5 Disinfection and sterilising

The instructions relating to cleaning, disinfection and sterilisation protocols for accessories supplied by SATELEC, a company of Acteon group have been approved for each medical device and accessory. The applicable guides are listed in the chapter *Associated documentation page 1*.

They can be downloaded at the following address: www.satelec.com/documents.

Download



Instructions For Use

In all cases, the local regulations in force relating to the cleaning, disinfection and sterilisation protocols for accessories take precedence over the information provided by SATELEC, a company of Acteon group.

6 Regulations and standards

6.1 Latest document update

12/2021

6.2 Manufacturer identification



SATELEC
A Company of the ACTEON Group
17, avenue Gustave Eiffel
ZI du Phare
33700 MERIGNAC
France
Tel. +33 (0) 556.34.06.07
Fax. +33 (0) 556.34.92.92
E-mail: satelec@acteongroup.com
www.acteongroup.com



6.3 Manufacturer responsibility

The manufacturer shall under no circumstances be liable in the following cases:

- Non-compliance with manufacturer recommendations during installation, whether this is the network voltage or the electromagnetic environment
- Maintenance or repair procedures performed by people who are unauthorised by the manufacturer.
- Use on an electrical fixture that is not compliant with regulations in force.
- Use of the device for purposes other than those specified in this manual.
- The use of accessories or handpieces other than those supplied by SATELEC, a company of Acteon group.
- Non-compliance with the instructions contained in this document.

| Note: the manufacturer reserves the right to modify the medical device and any documentation without notice.

6.4 Branch addresses

AUSTRALIA/NEW ZEALAND

ACTEON AUSTRALIA/NEW ZEALAND
Suite 119, 30-40 Harcourt Parade
Rosebery NSW 2018
Australia
Tel. +612 9669 477 307
Fax. +612/96692204
info.au@acteongroup.com

BRAZIL

MICRO IMAGEM INDUSTRIA COMERCIO IMPORTAÇÃO E EXPORTAÇÃO LTDA
CNPJ: 14.041.012/0001-79
Alameda Vênus, 233
Distrito Industrial
Indaiatuba – SP – CEP 13347-659
Brazil
Tel. +55 19 3936 809

CHINA

ACTEON CHINA
Office 401 - 12 Xinyuanxili Zhong Street -
Chaoyang District - BEIJING 100027 - CHINA
Tel. +86 10 646 570 11/2/3
Fax. +86 10 646 580 15
info.cn@acteongroup.com

RUSSIA

ACTEON RUSSIA
Gilyarovskogo str, 6b1, off 212
129090 MOSCOW
RUSSIA
Mob. +7 926 233 1695
Tel. +7 495 150 1323
info.ru@acteongroup.com

SPAIN

ACTEON MEDICO-DENTAL IBERICA, S.A.U.
Avda Principal n°11 H
Poligono Industrial Can Clapers
08181 SENTMENAT (BARCELONA) - SPAIN
Tel. +34 93 715 45 20
Fax. +34 93 715 32 29
info.es@acteongroup.com

TAIWAN

ACTEON TAIWAN
11F., No.1, Songzhi Rd.
Xinyi Dist., Taipei City 11047
TAIWAN (R.O.C.)
+ 886 2 8729 2103
info.tw@acteongroup.com

GERMANY

Acteon Germany GmbH
Klaus Bungert Straße 5, D-40468 Düsseldorf
GERMANY
T: +49 211 16 98 00-0
F: +49 211 16 98 00-48
www.acteongroup.com/de-de

INDIA

ACTEON INDIA
1202, PLOT NO. D-9
GOPAL HEIGHTS, NETAJI SUBASH PLACE
PITAMPURA, DELHI - 110034 - INDIA
Gujarat - India
Tel. +91 11 47 018 291 / 47 058 291 / 45 618 291
Fax. +91 79 2328 7480
info.in@acteongroup.com

ITALIA

ACTEON ITALIA
Via Roma 45
21057 OLGiate OLONA (VARESE)
ITALY
Tel. +39 0331 376 760
Fax. +39 0331 376 763
info.it@acteongroup.com

LATIN AMERICA

ACTEON LATINA AMERICA
Bogotá - COLOMBIA
Mobile: +57 312 377 8209
info.latam@acteongroup.com

MIDDLE EAST

ACTEON MIDDLE EAST
247 Wasfi Al Tal str.
401 AMMAN - JORDAN
Tel. +962 6 553 4401
Fax. +962 6 553 7833
info.me@acteongroup.com

THAILAND

ACTEON (THAILAND) LTD
23/45 Sorachai Building 16th floor - Sukumvit 63
Road, Klongton Nua - Wattana, BANGKOK 10110
- THAILAND
Tel. +66 2 714 3295
Fax. +66 2 714 3296
info.th@acteongroup.com

TURKEY

ACTEON TURKEY
Barbaros Mah. Dereboyu Cad.
Akzambak Sokak Uphill Towers
B Blok K14 D84 Atasehir
ISTANBUL 34746 - TURKEY
Mob. +90 532 481 20 57
Tel. +90 216 688 88 68
talha.gonca@acteongroup.com

U.K.

ACTEON UK
Phoenix Park– Eaton Socon, St Neots
CAMBS PE19 8EP - UK
Tel. +44 1480 477 307
Fax. +44 1480 477 381
info.uk@acteongroup.com

U.S.A. & Canada

ACTEON North America
124 Gaither Drive, Suite 140
Mount Laurel, NJ 08054 - USA
Tel. +1 856 222 9988
Fax. +1 856 222 4726
info.us@acteongroup.com

6.5 Disposal and recycling

As an item of Electrical and Electronic Equipment, the medical device must be disposed of via a specialist collection, removal, recycling or destruction channel. This applies in particular to the European market, with reference to Directive no. 2012/19/EC of July 2012.

When your medical device has reached the end of its service life, contact your nearest dental equipment dealer, or the Acteon head office or one of the company branches to find out how to proceed. The relevant contact details are given in the chapter *Branch addresses page 13*.



| The indication below applies to France only.

In compliance with the provisions of the French Environment Code relating to the disposal of electronic and electrical equipment waste or WEEE (Decree no. 2012-617 dated 2 May 2012), our Company fulfils its obligations to reclaim and dispose of its electrical and electronic equipment through the means established by the approved organisation Récyllum, NOR approval: DEVP1427651A.

As a manufacturer, our Company is listed in the National Register of Producers kept by the ADEME (French Environment and Energy Management Agency). Professionals buying our products directly from the distribution chain are responsible for passing on this information about our established recycling methods to the end user.

In addition, the buyer agrees to take back our brand's devices at the end of their service life and to transfer them to one of the collection centres set up by Récyllum for recycling (see list of collection centres on the site <http://www.recyllum.com/>).

If necessary, Récyllum can come and collect these devices from you free of charge once the quantity of devices has reached a certain level in the pallets-containers with which you are provided to store this waste.



A medical device that has reached the end of its service life must be disposed of in infectious clinical waste containers.

6.6 Applicable standards and regulations

This medical device complies with the essential requirements of European Directive 93/42/EEC. It was designed and manufactured in accordance with an EN ISO 13485-certified quality assurance system.

6.7 Symbols

| Symbol | Meaning |
|--|---------------------------------------|
|  | Always wear safety goggles |
|  | Always wear protective gloves |
|  Refer to Instruction Manual/Booklet | Refer to the supporting documentation |
|  Consult Instructions for Use | Consult the User Manual |

| Symbol | Meaning |
|---|--|
|  <p data-bbox="344 293 442 322">Electronic User Information</p> | The accompanying documentation is available in electronic format |
|  | Biohazard |
|  | Sterilisation at 134°C in an autoclave |
|  | Sterilisation at 132°C in an autoclave |
|  | Washer-disinfector for thermal disinfection |
|  | CE marking |
|  | Year of manufacture |
|  | Manufacturer |
|  <p data-bbox="339 1256 442 1308">Do not dispose of as household waste</p> | Do not dispose of as household waste |
|  | Recycle your lamps and professional electrical equipment with Réylum |
| <p data-bbox="336 1458 427 1487">Rx Only</p> | Under the United States Federal Law, this medical device must only be sold by or under the orders of a qualified doctor. |

7 Glossary

A

Active diameter

area of the effective optical cross-section of the LED light beam at the optical guide tip

alcohol wipe

disposable wipe soaked in an alcoholic solution designed to disinfect medical devices

autoclave

container with thick walls and hermetic seal designed to steam sterilise under a pressure of several bar. For an item to be considered sterile, the theoretical probability of isolating a germ must be less than 1 in a million. This is the sterility assurance level (SAL) stipulated in standard EN 556.

C

cleaning

essential pre-conditioning step to remove contamination through the physical-chemical action of a suitable product such as a detergent, combined with a mechanical action to ensure that the medical device is fully operational and clean. After cleaning, the cleanliness of the medical device components should be checked in addition to the cleanliness of the reassembled medical device. It is also important to make sure there is no damage likely to impact the safety, integrity or correct operation of the device

D

disinfection

voluntary and temporary removal of some germs to stop or prevent an infection, risk of infection or secondary infection by unwanted or pathogenic viruses or micro-organisms

E

expiry date

date up to which the medical device can be used. After this date, the medical device will need to be

resterilised

F

Ferrule

metal ring placed on the end of the optical guide. Makes it easier to insert the optical guide into the handpiece nozzle and prevents the optical guide from rotating.

Flexible protection shield

available in 5.5 mm-diameter and 7.5 mm-diameter sizes. In contact with the patient, it must be sterilised by autoclave before and after each use. Previously called the cup

I

Irradiance

term used in radiometry to quantify the power of an electromagnetic radiation per unit area. It is expressed in watts per square metre. Often confused with the power of a light source

L

LED

electroluminescent diode, more commonly known as Led (light-emitting diode). Designates an optoelectronic component that allows the emission of monochromatic light

O

operator

practitioner using the medical device during a treatment

Optical guide

light conductor fitted to the handpiece nosepiece and transmitting light to the cure site. Is cleaned, disinfected and sterilised in an autoclave.

P

practitioner

medically qualified person responsible for buying and operating the medical device

pre-disinfection

initial treatment to be performed on contaminated objects and equipment in order to reduce the number of micro-organisms and to facilitate subsequent cleaning. It is important to prevent residue from drying on the equipment. The other purpose of pre-disinfection is to protect personnel during the handling of instruments and to protect the environment. It is performed as soon as possible after use of the medical device within the vicinity of the place of use, prior to cleaning and in accordance with a procedure validated by the quality assurance system manager. The bactericidal, fungicidal and virucidal activities of the products used are determined in accordance with standards in force. These products are compatible with the medical devices to be handled and do not contain any substance known to be able to bind proteins

pre-vacuum

forced extraction of air from inside the autoclave sterilisation chamber

Protection plug

two plastic plugs used to protect the handpiece connectors and electronics during cleaning. One fits to the nosepiece and the other fits to the handpiece's electrical connectors

R

Rigid protection shield

removable oval shield forming an integral part of the handpiece once in place. Is cleaned with wipes. Not suitable for autoclaving

S

sterilisation

process used to kill potentially infectious viable or revivable germs in medicines or on medical devices. By definition, the sterility of a medical device is determined by a 1 in 1,000,000 probability of finding a viable or revivable germ on (or in) a product

U

ultrasonic tank

or ultrasonic cleaning. Rapid part cleaning or product dissolution procedure using the mechanical effect of ultrasonic waves

user

practitioner using the medical device to perform a clinical procedure. Also called operator

W

washer-disinfector

device designed to clean and disinfect batches of surgical instruments, anaesthetic accessories, earthenware, utensils, glassware and similar items. Generally works by washing with a detergent, thermally disinfecting and drying, sometimes by means of vacuum

Wavelength peak

maximum amplitude of a wavelength spectrum

8 Index

C

cataract 6
Chlorotetracycline 5
Class 2 safety goggles 6

D

disposal 15

E

electronic 1
electronic user instructions 1
Erythropoietic protoporphyria 5
European directive 15

L

lens 6

M

Methoxsalen 5

R

recycling 15
Récyllum 15
rest time 6
retina 6

S

solar urticaria 5

U

update 13
User Manual 1

W

warning 3



General Instructions | MINILED | J05102EN | V7 | (09) | 12/2021 | RG41EN040G

SATELEC | A Company of the ACTEON Group
17 av. Gustave Eiffel | ZI du Phare
33700 MERIGNAC | FRANCE
Tel. +33 (0) 556 34 06 07 | Fax. +33 (0) 556 34 92 92
satelec@acteongroup.com | www.acteongroup.com

