

# Completes the team for the Hygoclave 90.

Perfect accessories for smooth day-to-day practice.



# Everything that you need



## Hygoprint Plus

The automatic label printer for sterile product packaging enables the direct printing of program-specific labels.

- Ethernet connection allows flexible positioning
- Label size in accordance with guidelines - can be applied outside the sealing seam
- Bar code for scanning in digital data management
- Label specifications: device (REF and serial number), program, LOT date, operator (approval by), expiration date

Order number: 6046100035



## Labels

Self-adhesive labels for refilling for the Hygoprint Plus. Ideal for digital and analogue laboratory organisation.

- Incl. wax colour ribbons
- Packaging size: 2 rolls each consisting of 2,500 labels + 2 wax colour ribbons

Order number: 6046100130



## Hygotest

The Helix test body with 250 test strips for a regular, quick and simple Helix test.

- In accordance with EN 867-5 and ISO 11140-1
- PCD (Process Challenge Device) and 250 indicators
- Can be used as a batch control system and Bowie-Dick test
- Starter set with 50 indicators included in the initial autoclave delivery

Order number: 6050-600-03



### Trays

For practical transport and neat sorting of the equipment.

- In standard tray sizes
  - Made of highly corrosion-resistant stainless steel
  - Can be thermally disinfected
- Order number: 6046100112



### Software Tyscor Pulse

The software enables convenient protocol management as well as the direct display and saving of status messages.

- Network monitoring of the Hygiclave 90
  - Graphical display of the process steps inc. remaining running time
- Order number: 0949-150-00



# Convenient water treatment. For a perfectly flowing hygiene system.

With the Hygodem 90 and Hygopure 90.

## **The Hygodem 90. Compact reliability everywhere.**

- Ion exchange unit
- Convenient thanks to 3 ion exchanger cartridges
- Relief valve makes for simple maintenance
- Quick connector makes for simple installation
- Low space requirements
- Affordable purchase price



## **The Hygopure 90. High capacity for high demand.**

- Reverse osmosis system incl. water tap and storage tank
- Multiple filter stages + high-performance diaphragm with a service valve for checking its effectiveness
- Quick connector makes for simple installation
- High capacity
- Several litres of DI water can be removed at once
- Eco module allows water saving of up to 70%

