

You cannot buy aspirations or standards. But quality – that is a different matter.



Without exception, Orotol® plus and MD 555 cleaner repeatedly passed all material tests in 2017.* This confirms the consistently high quality of the market leader – made in Germany.

An unbeatable duo that helps suction systems and amalgam separators to retain their value:

Daily disinfection of the suction unit is a matter of routine in every surgery. Nonetheless, around half of all suction units suffer from reduced performance due to deposits and blockages. The reason: although most suction units are disinfected, no additional cleaning is performed. However, this aspect is crucial, after all, we are dealing with lots of aggressive substances here, such as blood, secretions, tooth material, filling materials or prophylaxis powders. These can form a tough sludge. This leads to blockages that can lead to a drop in performance and, in the worst case, a complete standstill of the surgery as a result of failure of the suction system. Not only can Orotol® plus and MD 555 cleaner prevent this, but at the same time they also offer extremely good compatibility with all the materials*. The two concentrates complement each other perfectly: while the alkaline disinfectant Orotol® plus rigorously tackles bacteria, fungi and viruses, the acidic MD 555 cleaner reliably removes deposits that are difficult to dissolve from all parts of the suction system, including from the pipes.

Our quality promise:


Orotol® plus and MD 555 cleaner

- Effective removal of biofilm and deposits from all components of the suction system and permanent prevention of their buildup if used on a regular basis
- Extremely gentle on materials (proved in regular material tests)
- Ensure that suction systems and amalgam separators retain their value as best possible
- Recommended by leading equipment manufacturers
- Market leaders in Europe – made in Germany


Orotol® plus reliably kills and inactivates microbes throughout the entire suction system (wide scope of action: bacteria incl. Tb, fungi and viruses incl. HBV, HCV, HIV).

MD 555 cleaner

reliably dissolves and removes incrustations of difficult-to-dissolve salts, lime and prophylaxis powders (including pearl products).



*Only when Orotol® plus and MD 555 cleaner were used together as a system was it possible to demonstrate material compatibility on all parts of the suction system during the annual tests.



Amalgam separator: Clear differences are also evident in the cassettes of the Dürr Dental amalgam separators.

For more information please visit www.duerrdental.com



The killer duo.

Hygiene power for extra long life of your suction unit.



Orotol® plus and MD 555 cleaner work together to protect all of the components in a suction unit against deposits, microbial contamination and encrustation. In addition, the new formulation of MD 555 cleaner is even more effective at dissolving even the most difficult-to-remove residues from prophylaxis pearl products. Regular use of both products will extend the service life of your suction unit.

For more information please visit www.duerrdental.com

**DÜRR
DENTAL**
THE BEST, BY DESIGN