

Materials



Paro Steribox silver

Sterilisation box with perforated base for sterilising and storing the Vector Paro handpiece and tool kits

2031-330-00



Scaler Steribox blue

Sterilisation box with perforated base for sterilising and storing the Vector Scaler handpiece and tool kits

2032-330-00



Vector Fluid polish

Polishing fluid with hydroxylapatite particles for smoothing and polishing; contents: 200 ml

CWZ510C2350



Vector/RinsEndo disinfectant

Ready-to-use disinfectant solution for cleaning and disinfecting hose connections, DGHM-listed; packaging size: 2.5 l bottle

CDZ501C6101



Vector cleaner

Cleaner for the Vector hose system
packaging size: 2.5 l bottle

CCA531A6150

Therapeutic concepts and instrumentation

The Vector system offers comprehensive possibilities for the treatment of periodontitis and periimplantitis as well as for professional teeth cleaning.

P007-625-02/DD-dd.de/1.0/01/101 Subject to change without notice

DÜRR DENTAL SE
Höfigheimer Str. 17
74321 Bietigheim-Bissingen
Germany
www.duerrendental.com
info@duerrdental.com

 **DÜRR
DENTAL**
THE BEST, BY DESIGN

Vector – instruments, materials and accessories



COMPRESSED AIR
SUCTION
IMAGING
DENTAL CARE
HYGIENE

 **DÜRR
DENTAL**
THE BEST, BY DESIGN

Periodontology and initial treatment



Tool kit Paro

Instrument case for 5 Paro instruments
2031-450-00



Paro Probe straight

Periodontal probe-equivalent instrument for the treatment of subgingival and supragingival labial, buccal, and oral surfaces
2030-151-01E Package contents: 3 pieces



Paro Lancet

Lancet-shaped instrument for the treatment of subgingival labial, buccal, and oral surfaces in pockets of up to 5 mm
2030-151-02E Package contents: 3 pieces



Paro Probe curved

Naber's probe-equivalent instrument for the treatment of interradicular spaces
2030-151-03E Package contents: 3 pieces



Paro Curette

Angled instrument for the treatment of subgingival and supragingival approximal surfaces
2030-151-04E Package contents: 3 pieces



Paro Probe Plus

Periodontal probe-equivalent instrument for the universal treatment of approximal, oral, and vestibular surfaces. Access depth of up to 5 mm in approximal area
2031-400-06E Package contents: 3 pieces

Periodontology, preservation, and periimplantitis treatment



Tool kit Recall/Implant

Instrument case for 3 Recall instruments
2031-460-00



Recall Probe straight CFK

Probe-shaped, semi-flexible instrument (carbon fibre material) for the treatment of subgingival and supragingival labial, buccal, and oral surfaces
2030-153-02E Package contents: 3 pieces




Recall Curette CFK

Curved, semi-flexible instrument (carbon fibre material) for the treatment of approximal surfaces
2030-153-05E Package contents: 3 pieces




Supra Probe flexible

Probe-shaped, highly flexible instrument for the gentle cleaning and rinsing of implant and tooth surfaces
2030-152-01E Package contents: 3 pieces



Periimplant soft

Flame-shaped, partly flexible, universally applicable instrument (carbon fibre material) for the gentle cleaning of sensitive implant surfaces.
2031-474-01E Package contents: 3 pieces



Periimplant hard

Conical-shaped instrument (carbon fibre material) for the reliable and efficient removal of dental plaque on implant surfaces.
2031-473-01E Package contents: 3 pieces

Prophylaxis



P4 Scaler instrument

Universal instrument for the removal of supragingival plaque
P4: 2032-414-00



P1, P2, P3 Scaler instruments

For the removal of subgingival plaque in deep gum pockets
P1: 2032-411-00 (straight)
P2: 2032-412-00 (curved to the right)
P3: 2032-413-00 (curved to the left)

Convince your patients!
We can provide free support



Vector marketing package

- Patient brochures with display
- Patient poster
- Prophylaxis cards

Vector DVD

- An application film explains how the instruments are used
- More documents and information about the Vector system