

#whdentalwerk



video.wh.com

Now at
your dealer
or **wh.com**



W&H couplings

Versatile system interfaces



rotor quick

Versatile system interfaces

The Roto Quick coupling is a connector that transfers the distribution of air, water, electricity and light between the supply hose of the dental unit and air drive.



Versatile
For all turbine handpieces, air motors and scalers with Roto Quick connections.



Compatible
Compatible with all standardized turbine hoses.



Hygienic
Non-retraction valve reduces the backflow of contaminated spray water and minimizes the risk of infection.



Intelligent
The power transmission enables you to take advantage of revolutionary technologies, such as the 5x Ring LED+.



Flexible
360° instrument rotation for easy and ergonomic handling.



Ergonomic
Lightest system available to date.



Customizable
Continuous regulation of the spray water supply for rapid adjustment by the user – even during treatment.



W&H Roto Quick couplings

Light, flexible and sterilizable:

W&H Roto Quick couplings are versatile and suitable for use with transmission instruments with LED+, LED and halogen lights. 24 months warranty.



Roto Quick couplings:	RQ-24	RQ-34	RQ-14	RQ-03	RQ-04
Light:	Yes	Yes	–	–	–
Spray regulator:	–	Yes	Yes	–	–
Connection:	6-hole	6-hole	Standard 4-hole	2(3)-hole	Standard 4-hole



W&H coupling for Multiflex® connection:

Light:	LED
Spray regulator:	–
Connection:	6-hole



W&H RM-34 LED coupling for Multiflex® connection

Simple and uncompromising: The W&H RM-34 LED coupling is ready for use right away. This means: The coupling with LED can be used in addition to conventional halogen couplings without making the slightest changes to the dental unit. Long lifespan – 24 months warranty.

Manufacturer:

W&H Dentalwerk Bürmoos GmbH
 Ignaz-Glaser-Straße 53, Postfach 1
 51111 Bürmoos, **Austria**
 t +43 6274 6236-0
 f +43 6274 6236-55
 office@wh.com
 wh.com

Distribution:

W&H (UK) Limited
 6 Stroud Wood Business Centre
 Park Street, St Albans, Hertfordshire AL2 2NJ, **United Kingdom**
 t +44 1727 874990
 f +44 1727 874628
 office.uk@wh.com
 wh.com

