



Mobile, powerful and stand-alone dental suction system DO M (Aspina) can fully replace the dental suction system on the dental unit. It provides high-performance, mobile, independent and ergonomic suction suitable even for demanding treatments. Suction hoses and holders with suction regulation are fully compatible with various suction tips. Separation and waste collection is realised into the built-in separation vessel. DO M is suitable for alternating work between multiple treatment rooms and is especially useful in places with no connection to the drainage system. For its quiet operation and high performance, it is popular everywhere, where high suction performance and the comfort of a patient are required.



DO M - with noise reduction cabinet

## Without dryer

	DO M
Dryer	✗
Noise reduction cabinet	✓
Suction	✓
Nominal voltage / frequency	230 V / 50 Hz 230 V / 60 Hz 110 V / 60 Hz
Duty cycle	S1 - 100 %
Operating voltage	24 V / 0.1 A DC/AC
Max. output capacity	1100 l.min <sup>-1</sup>
Max. vacuum	12 kPa
Noise level	≤ 48 dB (A)
Weight	36 kg

Dimensions

565 x 350 x 860 mm

- Data on weight and dimension are for informational purposes only
- Technical modifications are subject to change

