

# DK50 B



Compressor with a small 4-liter air tank and higher working pressure up to 11.5 bar, suitable for specific applications with limited space.



DK50 B - without noise reduction cabinet


DK50 BS - with noise reduction cabinet

## Without dryer

	DK50 B	DK50 BS
Dryer	✗	✗
Noise reduction cabinet	✗	✓
Nominal voltage / frequency	230 V / 50 Hz 230 V / 60 Hz	230 V / 50 Hz 230 V / 60 Hz
Motor performance	0.55 kW	0.55 kW
Air tank capacity	4 l	4 l
Compressor efficiency (working pressure 9,5 - 11 bar)	50 l.min <sup>-1</sup> / 10 bar	50 l.min <sup>-1</sup> / 10 bar
Duty cycle	S3 - 50 %	S3 - 50 %
Noise level	≤ 65 dB (A)	≤ 45 dB (A)
Weight	35 kg	47 kg
Dimensions	290 x 430 x 490 mm	380 x 525 x 575 mm

- Data on weight and dimension are for informational purposes only
- Technical modifications are subject to change

# Diagrams



# Photogallery

