

# Product information filters

## Filters suction systems



Exhaust air bacteria filter for VS 250 S,  
V/V/S/VSA 300 S, Variosuc and PTS 120

**Order no. 7120-143-00 (with mounting  
bracket and accessories)  
7119100010 (single filter)**



Exhaust air bacteria filter for V/V/S 600 - V/V/S 1200 S,  
V 2400, PTS 200 and clinic systems V 6000,  
V 9000, V 12000, V 15000 und V 18000

**Order no. 0705-991-50 (with accessories)  
0705-991-05 (single filter)**

For hygienic reason we recommend installing a bacteria filter in the exhaust air line.

The Università di Camerino has tested the exhaust bacteria filters with the test virus PhiX174 and confirmed a retention of greater than 99.99 %. The test virus PhiX174 has a size of 25–30 nm.

Filter quality HEPA-H13 according  
EN 1822-1:2009

## Filters compressed air



Sterile filter for contamination reduction  
of the compressed air

**Order no. 1640-981-00**



Bacteria filter for contamination reduction  
of the compressed air

**Order no. 1650100172\***

Dürr Dental sterile filters and bacteria filters retain 99.99 % of particles greater than 0.22 µm.

The Università di Camerino has tested the sterile filter with the test virus PhiX174 and confirmed a retention of greater than 99.99 %. The test virus PhiX174 has a size of 25–30 nm.

Filter quality HEPA-H13 according  
EN 1822-1:2009

\* Replaces the sterile filter order no. 1640-981-00. When ordering the bacteria filter order no. 1650100172, the sterile filter 1640-981-00 will be delivered until approx. Q4/2020.

# Filter overview compressed air

Compressors with drying unit						
Compressor	Pcs.	Pcs.	Pcs.	Pcs.		
Primo	1x		1x	1x	or	1x
Duo	2x		1x	1x	or	1x
Trio	3x		1x	1x	or	1x
Quattro	4x		1x	1x	or	1x
Duo Tandem 1 Aggr.	2x		1x	1x	or	1x
Duo Tandem 2 Aggr.	4x		1x	1x	or	1x
Quattro Tandem 1 Aggr.	4x		1x	1x	or	1x
Quattro Tandem 2 Aggr.	8x		2x	2x	or	2x
Quattro P20	16x		4x	4x	or	4x
Tornado 1		1x	1x	1x	or	1x
Tornado 2/2+		1x	1x	1x	or	1x
Tornado 4		2x	1x	1x	or	1x
PTS 120		1x	1x	1x	or	1x
PTS 200 1 Aggr.	2x		1x	1x	or	1x
PTS 200 2 Aggr.	4x		1x	1x	or	1x



Silver Airline  
Compressor  
suction filter

0832-982-00




Tornado  
Compressor  
suction filter

5180-982-00



Sinter filter  
(only for membrane drying unit)

1650-101-00



Fine filter

1610-121-00



Old: Sterile filter  
1640-981-00

New: bacteria filter  
1650100172