

simpler

ENERGIZE YOUR PRACTICE!



simpler

JUST WHAT YOU NEED.




8W
power


Easy
to use


T.O.P.
Tissue
Optimized
Pulsing


7"
display
touch


Full
colour
screen


Android
platform


WiFi
Direct


Bluetooth


Light
handpiece


Pedal
wireless

SIMPLER is a diode laser for dentistry essential for treatments in

- ENDODONTICS
- PERIODONTICS
- IMPLANT
- THERAPY
- WHITENING
- DESENSITIZATION
- SURGERY FIBROTIC TISSUE
- SURGERY NORMAL TISSUE
- SURGERY GRANULATION TISSUE

The **TOP** Tissue Optimized Pulse mode is a biotype specific configuration that modulates laser emission for perfect bloodless cutting on every soft tissue type, with a laser power that can reach up to 8W. **SIMPLER** is the affordable and indispensable instrument for every modern dental practice!

WITH SIMPLER IT'S EASIER.

Simpler is extraordinary for the treatment of many different conditions in the oral cavity: on soft tissue, hard tissue, in periodontics, endodontics and surgery. Laser therapy can be an opportunity for photobiomodulation, reducing dentinal sensitivity, and healing herpes and aphthae. In aesthetics it is very effective for tooth whitening.

The **Simpler** diode laser by **Doctor Smile** is the perfect choice for the dental practice that needs a professional, easy, precise and versatile instrument.

Simpler is perfect for making your work easier and at the same time giving your patients the minimally invasive treatments they desire. Fast and intuitive, with a full colour touch screen display and Android platform, you'll quickly discover how easy a daily routine of painless treatment can be.

Ready for a new experience?



simpler

**DOCTOR
SMILE**
INNOVATION HAS A NAME

ESSENTIAL EVERYDAY.

Simpler will become a part of your studio and will satisfy your numerous therapeutical needs. You'll realize quickly how indispensable and productive its presence has become.

Your patients will be satisfied with the great sensation that laser therapy gives, since you will use less anaesthetics and sessions will be faster and without vibration.

Simpler simply effective, efficient and safe.



Made in Italy



Simpler is a CE medical device entirely made in Italy.

Model	Doctor Smile Simpler
Code	LA9D0 003.3
Wavelength	980 nm
Max power	8W
Pulse mode	Continuous, pulsed, super-pulsed
Frequency	Up to 25 kHz
Size (WxHxD)	225x154x137 mm
Weight	1,3 kg
Medical /Laser class	II B / 4

Attestato con Sistema di Gestione per la Qualità certificato secondo



All LAMBDA S.p.A. medical laser are certified with CE trademark in accordance with EEC Directives for medical device: EEC Directives 93/42

A M DP 1170

Dealer



MORE INFO

• doctor-smile.com • info@doctor-smile.com

**DOCTOR SMILE IS
A TRADEMARK OF
LAMBDA SPA**

Via dell'Impresa, 1
36040 BRENDOLA (VI) Italy
T +39 0444.349165
F +39 0444.349954
info@lambdaspa.com
lambdaspa.com