

Documentation and approval

Aligned solutions for the documentation of the entire decontamination process





MELAtrace® Pro and MELAtrace®

Software solutions for documentation, approval and traceability

MELA*trace*® Pro and MELA*trace*® comply with the requirements for complete documentation in an easy, safe and paper-less way. The intuitive MELAG software solutions connect the steps of cleaning and disinfection with sterilization in a breeze.

Just a few clicks are required to provide the batches with approval and to archive the program logs as well as all approval decisions in a **tamper-proof decontamination report**.

An additional MELA*print*® 60 or MELA*print*® 80 generates bar code labels in order to ensure a **complete traceability**. After the application of an instrument, you are able to transfer batch information directly from the bar code label into the patient records - irrespectively of your deployed patient management software.

Find out more about the advantages of the new MELA*trace*® Pro and MELA*trace*® documentation and approval solutions in our product comparison.





The MELA*trace*® Pro tutorials assist you in the documentation and approval process. Find more video tutorials on our multimedia channel:

www.melag.com/en/multimedia





MELAtrace®

convinces with unique product highlights at no extra cost:

- Intuitive operating interface:
 - Approval of the washer-disinfector, steam sterilizer and Careclave® in only a few clicks
- Legal certainty:

The program logs and approval decisions are saved in a tamper-proof decontamination report and archived in MELA*trace*®

Traceability:

Print bar code labels for marking the instruments by using the optional MELA*print*® 60 or MELA*print*® 80

Push notifications:

Automated notifications about important information and software updates

MELAtrace® Pro

combines all the product highlights of MELAtrace® with a number of further innovative features:

Load management:

Document the load elements and standard batches for maximum transparency

Video tutorials:

Aligned video tutorials provide all the information required for the proper use of MELA*trace*® Pro

Preferential support:

Rely on a preferential remote support by our MELAG experts

Sealing logs:

Assign sealing logs with all the relevant process parameters to the sterilization batch

① Complete documentation:

In addition to the load and approval information, MELA*trace*® Pro archives the text and graphic logs in the decontamination report

Attractive license model for 12, 36 or 60 months

Test now: 3-month free trial

See for yourself: Download MELAtrace with a 3-month trial license from our website and break the boundaries of the traditional documentation process. Afterwards, benefit from an attractive license model for 12, 36 or 60 months and access to the latest software and feature updates!


Documentation and approvalThe ideal software for the MELAG System Solution

First class in both price and performance: Our software solutions cover the whole process of instrument decontamination to ensure safe documentation. Discover the perfect interplay of MELAtrace® Pro and MELAtrace® with the MELAG System Solution.



Our digital tools

Smart workflows in instrument decontamination



- Media Center: Access a range of helpful video tutorials and user manuals.
- Troubleshooting: Identify malfunctions and receive solutions on the smartphone and tablet.
- * MELA*connect* App can be run on iOS and Android smartphones and tablets. The app can be connected with MELA*therm*® 10 Evolution, the steam sterilizers of Premium-Class Evolution and Cliniclave® 45 series via network integration.



Download Center:

Retrieve MELAG documents from anywhere and at any time:





Multimedia:

Learn from more than 300 videos about safe instrument decontamination:



www.melag.com/en/ multimedia



Troubleshooting:

Enter malfunction numbers to access immediate solutions:



www.melag.com/en/service/ troubleshooting

Perfection in the finest detail

Solutions for checking, logging and labelling

Our accessories aligned to the MELAG System Solution support you in checking the process success, saving logs and labelling wrapped instruments. Use our versatile range of accessories to customize the workflow to your needs and to ensure the best results in infection control.



Routine and batch control



MELAcontrol® Helix

is a helix test body with 250 indicator strips which simulates the most difficult requirements placed on the reprocessing of hollow body instruments in a class B steam sterilizer.



MELAcontrol® Pro

is designed for largescale sterilization needs. The re-usable test body guarantees even more economical batch control.



MELAcontrol® Seal Check

verifies the quality of the seal seam: Faulty fusion of the sterilization pouch or insufficient contact pressure is revealed through the indicator paper.



MELAcontrol® Ink Test

is an ink test with 30 single-use bags for the reliable control of sterilization foils with folds for irregularities, defects or channels.



Log transfer and printout



MELA*flash*® CF card printer

easily saves logs on a CF card without connection of the MELAG devices to the clinic network.



MELAnet® Box

automatically saves logs on your computer or server using a network connection.



MELAprint® 44

prints logs in paper form. An additional network adapter enables connection with new device generations.



MELAG USB Stick

for the documentation of Careclave®, the latest autoclave generations and MELAseal® 200 log files.



Label printing



MELAprint® 60 / 80

prints bar code labels to ensure optimal instrument marking and traceability.



Bar code hand scanner

transfers information from bar code labels into the digital patient record software.



MELAdoc®

generates labels with all the important information about the sterilization process for marking the wrapped instruments.



Documentation sheets

facilitates approval of the washer-disinfector and steam sterilizer batch in paper form.

Rev.3-23/0538-02.23-EN-ME

Overview of compatibility

			- 0			- SSS	- 333 				
	MELA <i>therm</i> ® 10 Serie	Careclave®	Cliniclave® 45 Serie	Vacuclave® 550	Premium- Plus-Klasse	Pro Line	Prime Line	Euroklav® und MELA <i>quick</i> ® 12	MELA <i>seal</i> ® Pro	MELA <i>seal</i> ® 200	MELA <i>seal</i> ® 100+
Routine and batch control											
MELAcontrol® Helix Order no. ME01080		•	•	•	•	•	•				
MELAcontrol® Pro Order no. ME01075		•	•	•	•	•	•				
MELAcontrol® Seal Check Order no. ME01079									•	•	•
MELAcontrol® Ink Test Order no. ME01089									•	•	•
Log transfer and storage											
MELAflash® CF-Card-Printer Order no. ME01039	Function integrated		Function integrated		Function integrated			•	•		
MELAflash® Kartenlesegerät Order no. ME01048	•		•		•			•	•		
MELAflash® CF-Card Order no. ME01043	•		•		•			₽ ¹	₽ ¹		
MELAnet® Box Order no. ME40296	Function integrated	Function integrated	Function integrated	Function integrated	Function integrated	Function integrated	Function integrated	•			
USB/Seriell- Adapter Order no. ME80270								•	•		
USB-Stick Order no. ME19901		•		•		•	•			•	
Log printout											
MELAprint® 44 Order no. ME01144	✓ ²		2		✓ ²			•	•		
MELAG Netzwerkadapter Order no. ME40295	•		•		•						
Documentation and approval software											
MELA <i>trace</i> ®	•	•	•	•	•	✓ ³	✓ ³	•	•	•	
MELAtrace® Pro	•	•	•	•	•	✓ ³	✓ ³	•	•	•	
Barcode- Handscanner ⁴ Order no. ME01110	✓ ⁴	✓ ⁴	✓ ⁴	✓ ⁴	✓ ⁴	✓ ⁴	✓ ⁴	✓ ⁴			
Label printing											
MELAprint® 60 Order no. ME01160	₹ 5	•	•	•	•	✓ ⁵	•	⊘ ⁵			
MELAprint® 80 Order no. ME01108	✓ 5	•	✓ 5	•	✓ 5	✓ 5	•	✓ ⁵			
MELAdoc® Order no. ME01095		•	•	•	•	•	•	•			

¹ In connection with MELA*flash®* CF card printer, ² In connection with MELAG Network Adapter, ³ Fully compatible from 01.07.2023, ⁴ In connection with MELA*trace®* Pro or MELA*trace®* and a computer with an USB interface, ⁵ In connection with MELA*trace®* Pro or MELA*trace®* and a computer with RJ45/Ethernet interface

