



System overview adaptor Assistina 301 plus

Straight and contra-angle handpieces

All straight and contra-angle handpieces with ISO connection

Contra-angles with removable heads



Adaptor REF 02679000

Sirona T1 Classic straight and contra-angle handpieces



Adaptor REF 02691000

FG chuck system (W&H, Sirona)



Adaptor REF 02693000

Entran contra-angle head



Adaptor REF 05204600

Turbines

Sirona T1, T2, T3



Adaptor REF 02692000

W&H turbines, air scaler and air motors with Roto Quick system



W&H Adaptor REF 02690400

Turbines or adaptors with 4-hole fix connection (with or without light)



Adaptor REF 04713200

Turbines with Multiflex^{®***} connection

Turbines with Borden 2(3)-hole fix connection



Adaptor REF 07014500

Base adaptor
REF 02685000



for 4-hole turbines



***) Multiflex[®] is a trademark of a third party not affiliated with W&H Dentalwerk Bürmoos GmbH.