

The best team for two operators – Tyscor VS 2, Tornado 2+ and CA 2

network
compatible
with VistaSoft
Monitor



COMPRESSED AIR
SUCTION
IMAGING
DENTAL CARE
HYGIENE

Perfectly matched, the team of Tyscor VS 2, Tornado 2+ and amalgam separator CA 2 convinces with excellent performance and top efficiency for two operators.

- Perfectly matched components for two operators
- All systems equipped with network interface
- Comfortable installation – plug 'n' play
- Monitoring via software VistaSoft Monitor
- State of the art technology and excellent design
- Super Silent – Noise level < 60 dB(A)

Powerful, network-compatible, and perfectly matched

Compressor Tornado 2+

- Electronic control with monitoring function for highest operational reliability
- Plug 'n' play installation – place, plug in, ready
- Anytime flexibly upgradeable with 2nd Tornado 2+, in same compressed air net
- Highest reliability and extremely low-noise
- Delivery at 5 bar: 125 l/min



Suction Machine Tyscor VS 2

- Reliable – thanks to high-quality components
- Energy-efficient – energy saving of more than 50 %
- Compact and light
- Can be retrofitted to convert a dry machine into a wet machine
- An intelligent control allows operation with two Tyscor suction machines connected to only one suction line



Software VistaSoft Monitor

- Overview of operational status of all connected devices in the task bar of your PC
- Central management of the connected units
- Maintenance intervals are always under control, for the highest degree of operational safety
- Quick and safe diagnoses for service purposes

Amalgam separator CA 2

- Centrifugal separator with hydrodynamic self-cleaning function
- Separation rate of > 95 %
- Electronic fill level monitoring
- Pressure expansion tank and powder-coated console are integrated



DÜRR DENTAL SE
Höpfigheimer Str. 17
74321 Bietigheim-Bissingen
Germany
www.duerrdental.com
info@duerrdental.com

 **DÜRR
DENTAL**
THE BEST, BY DESIGN