



# VALLO<sup>®</sup>

C O R D L E S S

HAVE YOU SEEN  
**THE LIGHT?**

[www.valo-led.com](http://www.valo-led.com)

**ULTRADENT**  
PRODUCTS, INC.

*Improving Oral Health Globally*

# VALO<sup>®</sup>

C O R D L E S S



A revolutionary broadband LED curing light,  
**VALO Cordless** offers:

- An optimally collimated beam that delivers consistent, uniform power over a range of surfaces and working distances.
- Three curing modes—Standard, High, and Xtra Power—for maximum versatility.
- Unibody construction that yields a durable, slim, ergonomic shape and allows unprecedented access to all restoration sites.
- Highly efficient LEDs that keep the wand body cool to the touch.
- A battery-operated, cordless wand for optimal convenience and flexibility.
- Operates on environmentally responsible, safe, inexpensive, rechargeable batteries.

## ADDITIONAL INFORMATION:

### #5941 – VALO Cordless Kit – \$1,534.49

- 1 - VALO Cordless LED curing light
- 4 - Batteries (rechargeable)
- 1 - Battery charger
- 1 - Power supply with universal plugs for battery charger
- 1 - VALO Cordless light shield
- 1 - Handpiece bracket holder
- 50 - Barrier sleeves

### Accessories:

- #5962 - Battery Charger - \$29.99
- #5963 - Rechargeable Batteries (2 qty) - \$10.99
- #5961 - Power Supply with Universal Plugs for Battery Charger - \$29.99
- #5964 - Barrier Sleeve Refill (500 qty) - \$30.99
- #5934 - PointCure Lens\* (2 qty) - \$30.99
- #5936 - ProxiCure Ball Lens\* (2 qty) - \$30.99
- #5937 - TransLume Lens\* (2 qty: green and orange) - \$30.99
- #5938 - EndoGuide Lens\* (2 qty) - \$30.99
- #5939 - Black Light Lens\* (1 qty) - \$40.99
- #1667 - Surface Mounting Bracket - \$5.99
- #5929 - VALO Cordless Light Shield\* (1 qty)\* - \$19.49

\*Reusable. Should be disinfected using an intermediate-level disinfectant.

Learn more and read the testimonials at  
[valo-led.com](http://valo-led.com)

800.552.5512 [ULTRADENT.COM](http://ULTRADENT.COM)

© 2013 Ultradent Products, Inc. All Rights Reserved.  
88942.5 011113

**ULTRADENT**  
PRODUCTS, INC.

Improving Oral Health Globally